

Model SV005 Sieve Shaker

- Uses electromagnetic impulses to impart a triple vibrating action (vertical, lateral and rotational)
- Simple and sturdy construction making it ideal for continuous and intense use
- Operated and controlled via a separate control panel enabling high precision, performance and repeatability
- Will accept 200, 250, 300 and 315mm, 8 and 12 inch diameter sieves
- Can hold up to 10 sieves
- Digital 1 to 999 minute timer
- Digital control panel gives control of the timer, vibration intensity and pauses between sieving actions

Code	Dims, w x d x h, mm	Weight, kg	Price
SIE1202	440 x 380 x 1080	65	£1791.00



Woven Mesh Sieves

- Stainless steel woven wire mesh in seamless stainless steel frames
- Precise fit, added strength and superior quality
- Guaranteed to meet ISO 3310 and BS 410 standards
- Dims, Ø x d: 200 x 50mm

Code	Mesh, mm	Price
SIE1000	16.00	£84.00
SIE1002	13.20	£84.00
SIE1004	11.20	£84.00
SIE1006	9.50	£84.00
SIE1008	8.00	£84.00
SIE1010	6.70	£84.00
SIE1012	5.60	£84.00
SIE1014	4.75	£84.00
SIE1016	4.00	£84.00
SIE1018	3.35	£67.00
SIE1020	2.80	£67.00
SIE1022	2.36	£67.00
SIE1024	2.00	£67.00
SIE1026	1.70	£67.00
SIE1028	1.40	£67.00
SIE1030	1.18	£67.00
SIE1032	1.00	£67.00
SIE1034	0.850	£67.00
SIE1036	0.710	£67.00
SIE1038	0.600	£67.00
SIE1040	0.500	£67.00
SIE1042	0.425	£67.00
SIE1044	0.355	£67.00
SIE1046	0.300	£67.00
SIE1048	0.250	£67.00
SIE1050	0.212	£67.00
SIE1052	0.180	£67.00
SIE1054	0.150	£67.00
SIE1056	0.125	£67.00
SIE1058	0.106	£67.00
SIE1060	0.090	£67.00
SIE1062	0.075	£67.00
SIE1064	0.063	£67.00
SIE1066	0.053	£86.00
SIE1068	0.045	£119.00
SIE1070	0.038	£124.00

