


Corning® Disposable Anti Static Spatulas



- Designed to save time and to provide contamination free samples
- Eliminates the cleaning and re-sterilising necessary with reusable spatulas
- Individually packaged, sterile, RNase/DNase free
- Anti static
- Non pyrogenic
- Doubled ended, available in different configurations




Code	Description	End 1	End 2	Price
3003	Spatula	Tapered blade	Spoon	£244.00
3004	Spatula	Narrow spoon	Spoon	£244.00
3005	Spatula	Round end blade	Spoon	£244.00
3006	Spatula	V shaped scoop	Spoon	£244.00
3007	Spatula	Flat end blade	Spoon	£244.00
3012	Microspatula	V shaped scoop	Tapered blade	£136.00
3013	Microspatula	V shaped scoop	Rounded blade	£136.00

Heathrow Scientific Smart Spatulas™



- Size related designs to fit perfectly with standard 0.5, 1.5/2.0, 15 and 50mL tubes
- Ideal for dispensing or sample collection of powders or soft solids
- Shatter resistant in liquid nitrogen
- Resistant to dilute acids and bases
- Can be used with hot liquids
- Autoclavable
- Complies with FDA standards for food contact



Code	Alt Ref	Description	Sterility	For Tubes, mL	Dims, w x l, mm	Pack	Price
SPA3104	HS120005	Micro anti static	Non sterile	0.2 and 0.5	140 x 3.5	300	£40.70
SPA3105	HS120008	Eco	Non sterile	15 and 1.5	140 x 7	300	£56.00
SPA3100	HS120003	Standard blue	Non sterile	15 and 1.5	210 x 7	300	£57.00
SPA3102	HS120004	Standard opaque	Non sterile	15 and 1.5	210 x 7	300	£58.00
SPA3108	HS120007	Standard opaque	Sterile	15 and 1.5	210 x 7	100	£227.00
SPA3106	HS120006	Macro green	Non sterile	50	310 x 10	150	£66.00

Spoon Spatulas

Code	Overall Length, mm	Material	Price
SPA1110	120	Nickel	£11.00
SPA1112	150	Nickel	£13.00
SPA1114	180	Nickel	£16.00
SPA1111	120	Stainless steel	£4.00
SPA1113	150	Stainless steel	£5.00
SPA1115	180	Stainless steel	£7.00
SPA1118	190 (offset spoon)	Stainless steel	£10.41
SPA1117	210	Stainless steel	£11.00

