



Seward Stomacher® 400 Circulator

seward
World leaders in sample preparation

The Stomacher® 400 Circulator allows for easy and effective processing of multiple samples, without the risk of cross contamination. The unique patented paddle design improves organism recovery by enhancing the stirring and extrusion forces applied to the sample so driving more organisms into suspension for analysis.

- Improved organism recovery with the patented curved shape of the Circulator paddles over older square paddle designs
- No risk of cross contamination between samples
- Durable and reliable
- Control panel with programmable settings or 'hands free' auto run
- Countdown of processing time
- Removable door for easy cleaning of paddle chamber

Code	Alt Ref	Description	Price
MIX3011	0400/001/UK	Stomacher® 400 Circulator	£3195.00




SLS Select Blender Bags



- Ideal for use in any blender
- Double size (4mm) seal ensures all faults are totally visible, giving the added security of no leakage
- Tested in production at eight times the normal usage
- Gamma irradiated
- Standard bags are made from irradiated virgin food grade polythene
- Eco bags are manufactured from a new polymer which gives the bag the same strength as a standard bag, whilst being over 20% thinner
- Curved bags have a curved lower seal for Circulator® machines
- Strainer bags have a full width filter suspended inside with a mesh size guaranteed to allow the free movement of organisms through it, whilst stopping debris and fat through, making it easier to pipette the sample for testing

Code	Description	For use with	Thickness, µm	Dims, w x l, mm	Pack	Price
SLS3100	Standard	80mL models	65	100 x 150	2000	£79.00
SLS3102	Standard	400mL models	65	180 x 300	500	£72.00
SLS3104	Standard	3500mL models	70	380 x 500	250	£113.00
SLS3103	Curved	400mL circulator® models	65	180 x 300	500	£71.00
SLS3108	Strainer	400mL models	65	180 x 300	200	£112.00
SLS3109	Strainer curved	400mL circulator® models	65	180 x 300	200	£100.00
SLS4102	Eco	400mL models	50	180 x 300	2500	£208.00



SLS Select Pipette Straws



- Disposable pipette straws, specifically designed for use with the SLS Select range of blender bags
- Sterile
- Unbreakable and eliminate the hazard of broken glass
- Avoid exposure to infectious materials and contaminating samples

Supplied in packs of 2000

Code	Colour	Length, mm	Price
SLS3112	Blue tint	180	£98.00
SLS3114	Clear	190	£89.00