

Spex® Jenway® 72 Series Spectrophotometers Continued



Accessories

Code	Alt Ref	Description	Price
SPE2234	83056-79	External printer	£461.00
CON1236	37702	Replacement paper roll for printer	£8.00
SPE2232	725201	Heated cell 10 x 10 cuvette holder with UK, EU and US power leads	£522.00
SPE3514	630204	Single cuvette holder, 10 x 10mm	£147.00
SPE3518	630005	Adjustable path length cuvette holder, 10 to 100mm	£147.00
SPE3520	630304	Micro cuvette holder with adjustable aperture	£146.00
SPR4682	637071	Test tube holder, 16/24mm	£277.00

Spex® Jenway® 73 Series Spectrophotometers



These 73 series single beam spectrophotometers have improved optics resulting in a narrow spectral bandwidth of 5nm and absorbance range of -0.3 to 2.5A. All models feature measurement modes for absorbance, % transmittance, and concentration. Advanced models add spectrum scanning, quantitation, and kinetics.

- Icon driven software and easy to use navigation system provide intuitive usability
- Autologging capabilities enable repeat measurements of the same sample with a selected time interval
- Minimal footprint with large graphical display
- Bandwidth: 5nm
- Absorbance: -0.3 to 2.5A
- Wavelength accuracy: ± 2.0 nm
- Wavelength reproducibility: ± 0.5 nm
- Dims, w x d x h: 280 x 220 x 400mm



3 year warranty, 1 year on the tungsten lamp

Code	Alt Ref	Model	Description	Light Source	Wavelength Range, nm	Price
SPE4600	730001	7300	Visible scanning spectrophotometer	Tungsten halogen	320 - 1000	£2488.00
SPE4620	730501	7305	UV/visible scanning spectrophotometer	Xenon	198 - 1000	£3162.00
SPE4640	731001	7310	Advanced visible scanning spectrophotometer	Tungsten halogen	320 - 1000	£2870.00
SPE4660	731501	7315	Advanced UV/visible scanning spectrophotometer	Xenon	198 - 1000	£3817.00

Accessories

Code	Alt Ref	Description	Price
SPE4668	660101	Internal printer	£464.00
SPE4670	735401	Automatic 8 position cuvette holder	£1305.00
SPE4672	735201	Sipper pump	£1471.00
SPE4674	735301	Peltier	£905.00
SPE4676	735701	Combined sipper/peltier	£1732.00
SPE3514	630204	Single cuvette holder, 10 x 10mm	£147.00
SPE3518	630005	Adjustable path length cuvette holder, 10 to 100mm	£147.00
SPR4682	637071	Test tube holder, 16/24mm	£277.00
SPE3520	630304	Micro cuvette holder with adjustable aperture	£146.00
SPR4686	736201	Water heated 10 x 10mm path length cuvette holder	£525.00
SPR4688	735601	Boiling tube holder	£237.00
SPR4690	735001	Dust cover	£87.00