

Bruker BRAVO Handheld Raman Spectrometer Continued



Code	Alt Ref	Description	Price
SPE1076	BRAVO	BRAVO handheld Raman spectrometer	£35700.00

Bruker MultiRAM FT-Raman Spectrometer



The MultiRAM configuration provides a spectral range of 3600 – 50 cm⁻¹ (Stokes shift). The heart of the MultiRAM is Bruker's patented interferometer with gold coated optics and permanent RockSolid™ alignment, which provides the superb stability and throughput required for demanding experiments, as it is the case for Raman spectroscopy.

- Benchmark Raman shift accuracy and spectral resolution
- Large sample compartment for a variety of accessories
- High Throughput Sampling (HTS) and fibre coupling options available
- Microscopy coupling to RamanScopeIII as option
- Complies with GMP, 21 CFR part 11 & data integrity requirements

Prices are an approximate guide only. Contact us for a quote and to discuss configuration options.



Code	Alt Ref	Description	Price
SPE1078	MultiRam	MultiRAM FT-Raman spectrometer	£79700.00

Bruker TANGO Series FT-NIR Analysers



Faster, simpler, more secure – with TANGO your NIR analysis speeds up. TANGO has exactly what users require of an FT-NIR spectrometer suitable for industrial use: robustness, high precision and straightforward operator guidance. Bruker's proven technology was combined with an easy to use touch screen operation and a small footprint, perfect for those laboratories with limited space.

- Intuitive user interface available in 17 languages makes TANGO easy to understand
- Robust and precise optics
- Direct calibration transfer to and from all Bruker Optics FT-NIR spectrometer
- User friendly maintenance
- Data exchange via network
- Easy sample handling avoiding time consuming sample preparation, chemicals or additional reagents
- Complies with 21 CFR part 11 requirements

Prices are an approximate guide only. Contact us for a quote and to discuss configuration options.



Code	Alt Ref	Description	Price
SPE1080	TANGO-T	TANGO-T FT-NIR analyser for transmission measurements	£42300.00
SPE1082	TANGO-R	TANGO-R FT-NIR analyser for reflection measurements	£42300.00