

Seward Simplette® Straws

seward
World leaders in sample preparation

- Provide an accurate and economical alternative to pipette tips
- Designed specifically to pipette samples from the bottom of the Stomacher® bag without risk of cross contamination
- Capable of pipetting volumes between 0.1 – 2mL
- Single use
- Irradiated sterile



Code	Alt Ref	Length, mm	Inner Ø, mm	For use with	Pack	Price
SEW3112	BA6051	177	4.1 - 4.2	Barky µltipette™	2000	£74.00
SEW7114	BA6052	250	4.1 - 4.2	Barky µltipette™	1000	£102.00
SOC3112	BA6053	177	3.6 - 3.7	Socorex Acura 810	2000	£75.00

Barky SPR µltipettes™ Pipettes

- Specially designed pipettes for use with straw tips
- Plastic straw tips are a safe alternative to glass pipette where fast moving sampling or diluting is required
- Robust construction ensures a long life and reliable operation
- Stainless nozzle will take a straw tip securely avoiding fall offs
- Colour coded for easy identification in a busy lab
- Re-calibration if required at any time can be carried out in house
- Can be autoclaved as a unit or individual parts can be autoclaved
- Easy maintenance and servicing which can be carried out in house
- Smooth action of the plunger helps avoiding thumb fatigue
- Designed so hand warming effects are kept to the minimum
- Can be traced back through a serial number engraved on each instrument



Code	Description	Volume, mL	Price
PIP1500	SPR-1 variable volume	0.1 - 1.0 (increments of 0.1)	£185.00
PIP1510	SPR-100 single volume	0.1	£185.00
PIP1520	SPR-500 single volume	0.5	£175.00
PIP1530	SPR-1000 single volume	1.0	£175.00

View our Safety & PPE catalogue

To browse our new Safety and PPE catalogue, visit the SLS website:

scientificlabs.co.uk/literature (UK & Northern Ireland)

scientificlabs.ie/literature (Republic of Ireland)

scientificlabs.com/literature (International)

