





SLS Flowgen Mini Digital Dry Bath



With its compact construction, broad temperature range, and multiple block options, the mini digital dry bath is the perfect tool for laboratory incubations. The simple touch pad controls and a digital display provide for "set and walk away" temperature selection and unrivalled accuracy.

- Compact, fits in the palm of your hand (less than 12cm wide)
- Exchangeable blocks for 0.2 to 50mL tubes (sold separately)
- Clear cover ensures temperature uniformity
- Simple touchpad, digital control
- Supplied with block lifter

Code	Temp Range, °C	Accuracy, °C	Uniformity, °C	Increments, °C	Dims, w x d x h, mm	Price
BL02200	Ambient +5 to 100	±0.5	±0.2	0.1	105 x 150 x 115	£425.00



Blocks

High grade aluminium blocks, which are precision machined to match the conical size and shape of several different tube sizes. Assorted blocks are available for compatibility with many common laboratory tubes and plates, all of which provide unsurpassed temperature transfer.

Code	For use with	Price
BL02202	40 x 0.2mL tubes or 5 x PCR strips	£132.00
BL02204	24 x 0.5mL tubes	£101.00
BL02206	15 x 1.5mL tubes	£101.00
BL02208	15 x 1.5/2.0mL tubes	£101.00
BL02210	4 x 15mL tubes	£128.00
BL02212	2 x 50mL tubes	£128.00



SLS Flowgen Digital Dry Baths



Provide accurate "set and walk away" digital temperature selection, eliminating the need for external thermometers and repetitive fine tuning of a temperature control knob. Simply enter the desired temperature on the digital touchpad and the dry bath provides accurate temperature with real time monitoring that continuously maintains the selected temperature. A digital timer is built into the touchpad display, that can be used to alert the user that a chosen time interval has timed out.

- Include an easy to use block lifter with a temperature insulated handle
- Built in user re-calibration function, to guarantee long term temperature accuracy
- Temperature range: ambient +5 to 150°C
- Accuracy: ±0.2°C
- Uniformity: ±0.2°C
- Increments: 0.1°C

Code	No. Block Positions	Dims, w x d x h, mm	Price
BL02250	1	200 x 230 x 80	£526.00
BL02252	2	220 x 260 x 80	£656.00