


PTFE Covered Stirrer Magnets Continued



Pivot Ring

- For use in containers with bases that are slightly curved or uneven - the pivot enables the stirrer to adopt the optimum position for stirring

Code	Dims, Ø x h, mm	Price
STI5010	3 x 8	£0.80
STI5012	4.5 x 12	£0.83
STI5014	4.5 x 15	£0.83
STI5016	6 x 20	£0.96
STI5018	6 x 25	£1.17
STI5020	6 x 30	£1.17
STI5022	6 x 35	£1.29
STI5024	8 x 40	£1.74
STI5026	8 x 45	£1.69
STI5028	8 x 50	£2.09
STI5030	10 x 60	£2.33
STI5034	10 x 70	£3.00
STI5036	13 x 75	£4.98
IZ0013126	6 x 12	£0.90



Octahedral

- Similar action to pivot ring type

Code	Dims, Ø x h, mm	Price
STI4200	3 x 13	£1.28
STI4202	8 x 13	£0.73
STI4204	10 x 13	£1.10
STI4206	8 x 15	£1.00
STI4210	8 x 22	£1.24
STI4212	8 x 25	£1.24
STI4214	10 x 25	£1.42
STI4216	8 x 28	£1.49
STI4220	8 x 38	£2.00
STI4224	8 x 41	£2.15
STI4228	10 x 51	£2.26
STI4230	8 x 64	£3.02
STI4232	10 x 64	£2.82
STI4234	13 x 75	£4.00

PTFE Covered Rare Earth Stirrer Magnets

- PTFE covered stirrer bars with rare earth samarium-cobalt magnet cores
- Inert and leak free
- Isostatically moulded in pure, FDA approved PTFE
- Individually polished to stay cleaner for longer

- Exceptional strength
- Almost total resistance to demagnetisation
- Identified by an inert carbon black spot

See the SLS website for a larger range of shapes/sizes