




Shafts

- Inert, unbreakable, will not scratch
- Usable up to 280°C

Code	Shaft Ø, mm	Length, mm	Price
STI4300	6	400	£22.38
STI4302	6	500	£22.78
STI4304	9.5	550	£28.72
STI4306	9.5	650	£29.23





Blades

- Fit all standard stirring shafts
- Pure PTFE
- Totally inert and non scratch
- Hole diameter: 6.5mm
- Thickness: 3.2mm

Code	Description	Dims, w x h, mm	Price
STI4280	Square end	52 x 14	£1.85
STI4282	Square end	76 x 19	£2.40
STI4284	Square end	90 x 28	£3.28
STI4286	Plain end	65 x 25	£2.52
STI4288	Plain end	75 x 25	£2.68
STI4290	Plain end	105 x 25	£3.02
STI4292	Plain end	125 x 25	£3.19
STI4294	Plain end	150 x 25	£3.50

PTFE Stirrer Shafts

- Chemically resistant
- Inert
- Leak free
- Encapsulated stainless steel core for maximum rigidity
- Usable up to 280°C



Code	Description	Shaft Ø, mm	Length, mm	Rotor Ø, mm	Price
STI4330	Screw propeller	6	400	40	£39.78
STI4332	Screw propeller	6	500	50	£42.36
STI4334	Screw propeller	9.5	550	60	£60.20
STI4336	Screw propeller	9.5	650	60	£62.36
STI4340	Anchor propeller	6	400	80	£32.57
STI4346	Anchor propeller	9.5	550	140	£46.87
STI4348	Anchor propeller	9.5	650	140	£47.52
STI4352	Centrifugal shaft	6	500	50	£41.26
STI4354	Centrifugal shaft	9.5	550	70	£48.94
STI4356	Centrifugal shaft	9.5	650	70	£49.90