

SLS Select Pro bact Transport Swabs



- Complete M40 sampling kits for timed transport of all types of microorganisms
- CE marked
- Complies with the European Medical Device Directive 93/42/EEC
- Tamper evident seals and caps are included for assured security
- Sterilised by gamma irradiation to ensure a SAL of at least 10^{-6}



Amies Clear

- Gives up to 24 hours survival
- Ideal for microscopy work
- General clinical swabbing of fastidious, non fastidious and anaerobic microorganisms

Supplied in packs of 100

Code	Description	Price
SWA2062	Woodstick shaft with cotton tip	£70.00
SWA2064	Aluminium shaft with cotton mini tip	£110.00
SWA3074	Polystyrene (45mm breakpoint) shaft with viscose tip	£77.00
SWA2066	Polystyrene shaft with viscose tip	£64.00
SWA1110	Polystyrene shaft with cotton tip (3")	£105.00



Amies Charcoal

- Gives up to 48 hours survival using activated charcoal to absorb toxic metabolites away from the sample
- Ideal for clinically significant samples including fastidious and non fastidious microorganisms

Supplied in packs of 100

Code	Description	Price
SWA2068	Woodstick shaft with cotton tip	£70.00
SWA2071	Aluminium shaft with viscose tip	£111.00
SWA2072	Polystyrene shaft with viscose tip	£55.00
SWA3080	Polystyrene (45mm breakpoint) shaft with viscose tip	£77.00
SWA1100	Plasticised paper shaft with cotton mini tip	£102.00
SWA1104	Polystyrene shaft with cotton tip (3")	£105.00

SLS Select Enviroscreen Environmental and Hygiene Sampling Swabs



- A range of products to optimise hygiene and environmental monitoring protocols
- Sterilised by gamma irradiation to ensure a SAL of at least 10^{-6}



Pre-Moistened Swab Kits

- Convenient and simple sampling method to confirm the presence of microorganisms on equipment and critical surfaces
- No peel pouch, ready to use for convenience
- Pre-moistened to save time and cost
- Neutralise disinfectant residues
- Entire sample can be analysed
- Integral solution for HACCP
- Useful in food security/TACCP validation
- Supplied in tube

Code	Description	Pre-Moistened with	Pack	Price
SWA1002	Polystyrene (breakable) shaft with viscose tip	Neutralising buffer	250	£91.00
SWA1004	Polystyrene (45mm breakpoint) shaft with viscose tip	Neutralising buffer	250	£92.00
SWA1006	Polystyrene shaft with viscose tip	Deionised water	100	£33.00