



SLS Select Enviroscreen Environmental and Hygiene Sampling Swabs Continued



Dry Swabs

- Swabs for general sampling applications
- Supplied in a tamper evident sealed tube, with a write on label or in a peelable pouch

Code	Description	Packaging	Pack	Price
SWA2024	Polypropylene (flexible) shaft with foam tip (12")	Tube	100	£163.00
SWA1028	Woodstick with cotton tip	Tube	250	£82.00
SWA1030	Polypropylene (flexible) shaft with viscose tip	Tube	250	£88.00
SWA2028	Polystyrene (45mm breakpoint) shaft with viscose tip	Peel pouch	1000	£200.00
SWA1032	Polystyrene (45mm breakpoint) shaft with viscose tip	Tube	250	£87.00
SWA1034	Low ATP swab with a polystyrene shaft and polyester (dacron) tip	Tube	250	£109.00



Sampling Templates

- Sterile plastic template
- Incorporates press plate for easy sampling
- Can only be used once or crossover contamination will occur

Code	Dims, w x l, mm	Pack	Price
SWA3096	100 x 100	40	£75.00



3M™ Clean-Trace™ Hygiene Monitoring and Management System

3M Science.
Applied to Life.™

A simple, quick and reliable hand held diagnostic tool to verify cleaning effectiveness and capture, store and manage data for later analysis. It offers food safety professionals a complete system of products consisting of a luminometer handheld device, compatible tests and a powerful new software platform, to verify that surfaces have been effectively cleaned in seconds, for better decision making and audit preparedness. The system comes complete with luminometer and software and is used with the 3M™ Clean-Trace™ ATP tests to determine the level of contamination in a sample. The software turns data into information, helping you identify root causes and take action quickly.

- Re-engineered for improved precision and accuracy
- Easy to use and fast time to result
- Robust for the manufacturing environment
- Ergonomic design for simple one handed operation
- Wireless data transfer
- 23 languages for multilingual pre-production teams
- Robust to provide consistent results across temperatures
- Software has easy to use dashboard that immediately turns data into actionable information
- Includes the 3M™ Clean-Trace™ luminometer and 3M™ Clean-Trace™ hygiene monitoring and management software to help you monitor testing results over time (computer not included)

Code	Alt Ref	Price
FSA1260	LM1	£5389.00