



CLEARLock™ Microtubes

ClearLine

- Polypropylene microtubes combine the security of a screw cap tube with the ease-of-use of a snap cap tube.
- Designed to withstand up to 30,000 x g during centrifugation without any rupture and to preserve the sample integrity
- Hinged lid prevents any unintentional opening during use
- Optimal transparency
- Temperature range: -80 to +100°C
- Autoclavable (open) at +121°C for 20 minutes
- Graduated with writing area
- Manufactured in Clean Room Class ISO 8 (UNI EN ISO 14644-1) Class 100000 (US FED STD209E)
- Human DNA, DNase, RNase, Pyrogen, ATP, PCR inhibitors free certified
- CE-IVD marked

Code	Alt Ref	Volume, mL	Description	Inner Pack	Pack	Price
TUB1042	CL05.050.0001S	0.5	Sterile	1	50	£26.70
TUB1034	CL05.001.0500	0.5	Non sterile	500	500	£12.23
TUB1052	CL05.010.1000B	0.5	Non sterile bulk	1000	10000	£149.00
TUB1044	CL15.100.0001S	1.5	Sterile	1	100	£49.60
TUB1036	CL15.002.0500	1.5	Non sterile	500	1000	£25.40
TUB1054	CL15.010.1000B	1.5	Non sterile bulk	1000	10000	£184.00
TUB1046	CL20.100.0001S	2.0	Sterile	1	100	£49.90
TUB1038	CL20.002.0500	2.0	Non sterile	500	1000	£27.20
TUB1056	CL20.008.1000B	2.0	Non sterile bulk	1000	8000	£167.00
TUB1048	CL50.050.0001S	5.0	Sterile	1	50	£27.80
TUB1050	CL50.010.0020S	5.0	Sterile	20	200	£65.00
TUB1040	CL50.002.0100	5.0	Non sterile	100	200	£28.50
TUB1058	CL50.006.0500B	5.0	Non sterile bulk	500	3000	£276.00

CLEARLine® Screw Cap Microtubes

ClearLine

- Manufactured from medical-grade raw material to exceed the exacting standards demanded by researchers
- Three styles of screw cap to choose from and a choice of conical bottom or self-standing tubes
- Withstand high speed centrifugation up to 16400RPM or 25000xg and temperatures from -85 to +121°C
- Autoclavable to 121°C for 20 minutes
- Tubes made in medical polypropylene and cap in medical polyethylene
- Certified DNase/RNase/Pyrogen free to satisfy the highest demands in several biotechnological applications



Sterile Tubes with O-Ring Caps

Supplied in packs of 1000 (inner packs of 500)

Code	Alt Ref	Capacity, mL	Bottom	Height, mm	Collar	Writing Area	Price
TUB1060	CL604/S	0.5	Conical	47.1	No	No	£93.00
TUB1062	CL606/S	0.5	Conical	48.5	Yes	No	£119.00
TUB1064	CL614/S	0.5	Self standing	47.5	No	No	£93.00
TUB1066	CL616/S	0.5	Self standing	49.2	Yes	No	£119.00
TUB1068	CL724/S	1.5	Conical	47.1	No	Yes	£101.00
TUB1070	CL726/S	1.5	Conical	48.2	Yes	Yes	£117.00
TUB1072	CL734/S	1.5	Self standing	47.6	No	Yes	£101.00
TUB1074	CL736/S	1.5	Self standing	49.0	Yes	Yes	£117.00
TUB1076	CL744/S	2.0	Conical	47.2	No	Yes	£98.00
TUB1078	CL746/S	2.0	Conical	48.4	Yes	Yes	£125.00
TUB1080	CL754/S	2.0	Self standing	48.6	No	Yes	£98.00
TUB1082	CL756/S	2.0	Self standing	49.8	Yes	Yes	£125.00