



## Miele G7825 Washer Disinfector

**Miele**

A medium sized front loading washer-disinfector with bottom hinged door and a capacity of 10 DIN trays per batch.

- Single chamber system for cleaning, disinfection and drying
- Double insulation prevents heat and noise emissions
- Hygienic fresh water rinsing system with a change of water after every rinsing phase
- 2 spray arms and injector system for thorough cleaning
- 17 standard cleaning and disinfection programmes with 39 customised programme slots
- Electrical door locking
- Programme failure check
- Visual and acoustic signal at the end of the programme
- 2 sensors for temperature monitoring and control
- Dims, w x d x h, mm: 900 x 750 x 1974

Chamber Dims.	Hot Air Drying	Deion. Water Pump
541 x 617 x 683	Yes	Accessory



**Note:** Services (water, drainage and 3-phase power) need to be prepared prior to delivery, installation and commissioning.

Code	Model
WAS7000	G7825

### Accessories

Code	Description
WAS7002	E741/1 mobile injector unit
WAS7004	E742 module
WAS7006	E743 injector module
WAS7008	E744 injector module
WAS7010	E745/1 injector module
WAS7012	E746 injector module
WAS7014	E747 injector Module
WAS7016	E755 injector module
WAS7018	E757 mobile injector unit
WAS7020	E775 mobile unit
WAS7022	E735/1 mobile unit
WAS7024	E701 mobile unit
WAS7026	E702 top insert
WAS7028	A7/1 insert

## Just a click away...

Our website aims to improve our customer's online experience. With over 300,000 products available and easy navigation, visit today for more information.

[www.scientificlabs.co.uk](http://www.scientificlabs.co.uk)

