

Sartorius Glass Microfibre Filters Continued

SARTORIUS

Grade MGC

Glass microfibre filter grade for the analysis of suspended solids in wastewater according to EN 872 with a particle retention of 1.2µm, weight 54g/m² and thickness 0.24mm. Alternative grade to Whatman GF/C.

Supplied in packs of 100

Code	Alt Ref	Ø, mm	Price
SAR4304	FT-3-1103-025	25	£27.00
SAR4310	FT-3-1103-047	47	£40.00
SAR4311	FT-3-1103-050	50	£46.00
SAR4312	FT-3-1103-055	55	£47.00
SAR4314	FT-3-1103-070	70	£55.00
SAR4316	FT-3-1103-090	90	£80.00
SAR4317	FT-3-1103-100	100	£92.00
SAR4318	FT-3-1103-110	110	£105.00



Grade MGD

Prefilter with a particle retention of 2.7µm, weight 120g/m² and thickness 0.47mm. Alternative grade to Whatman GF/D.

Supplied in packs of 50

Code	Alt Ref	Ø, mm	Price
SAR4408	FT-3-1104-025	25	£15.00
SAR4414	FT-3-1104-047	47	£21.00
SAR4418	FT-3-1104-070	70	£34.00
SAR4420	FT-3-1104-090	90	£54.00
SAR4421	FT-3-1104-100	100	£63.00
SAR4424	FT-3-1104-125	125	£92.00



Our Life Science Brand

Equipment,
consumables
and reagents
for life science

