

## Whatman Puradisc Syringe Filters Continued

### 13mm Syringe Filters

- Sample volume up to 10mL
- Low hold-up volume < 25µL ensures maximum sample recovery
- Optional glass microfibre filter available
- Optional tube-tip format available
- Non-sterile

Code	Alt Ref	Membrane	Pore Size, µm	Pack	Price
FIL2740	67791302	PVDF	0.2	100	£180.00
FIL2742	67791304	PVDF	0.45	100	£180.00
FIL2744	67921302	PVDF	0.2	500	£655.00
FIL2746	67921304	PVDF	0.45	500	£644.00
FIL2748	67891301	Nylon	0.1	100	£177.00
FIL2750	67891302	Nylon	0.2	100	£177.00
FIL2754	67751302	PTFE with tube-tip	0.2	50	£120.00
FIL2756	67751304	PTFE with tube-tip	0.45	50	£117.00
FIL2758	67841301	PTFE	0.1	100	£177.00
FIL2760	67841302	PTFE	0.2	100	£180.00
FIL2764	67841310	PTFE	1.0	100	£177.00
FIL2766	67841350	PTFE	5.0	100	£172.00
FIL2768	67831302	PTFE	0.2	500	£649.00
FIL2770	67831304	PTFE	0.45	500	£644.00
FIL2772	67821302	PES	0.2	100	£180.00
FIL2774	67821304	PES	0.45	100	£180.00
FIL2776	67811304	PES	0.45	500	£758.00
FIL2778	67881302	PP	0.2	100	£177.00
FIL2780	67881304	PP	0.45	100	£177.00
FIL2782	67851302	PP	0.2	500	£649.00
FIL2784	67851304	PP	0.45	500	£649.00
FIL2786	67711304	CA	0.45	100	£172.00
FIL2788	68201316	GMF GF/A*	1.6	100	£177.00
FIL2790	68211310	GMF GF/B*	1.0	100	£177.00
FIL2792	68221312	GMF GF/C*	1.2	100	£177.00
FIL2794	68231327	GMF GF/D*	2.7	100	£177.00
FIL2796	68251307	GMF GF/F*	0.7	100	£177.00
FIL2798	68271315	GMF 934-AH*	1.5	100	£177.00




\* Particle retention rating

### 25mm Syringe Filters

- Sample volume up to 100mL

Code	Alt Ref	Sterility	Membrane	Pore Size, µm	Pack	Price
FIL5000	67802502	Yes	PES	0.2	50	£92.00
FIL5002	67812502	No	PES	0.2	200	£498.00
FIL5004	67802504	Yes	PES	0.45	50	£91.00
FIL5006	67812504	No	PES	0.45	200	£424.00
FIL5012	67842501	No	PTFE	0.1	50	£112.00
FIL5014	67842502	No	PTFE	0.2	50	£112.00
FIL5016	67852502	No	PTFE	0.2	200	£424.00
FIL5018	67842504	No	PTFE	0.45	50	£114.00
FIL5020	67852504	No	PTFE	0.45	200	£424.00
FIL5022	67842510	No	PTFE	1.0	50	£112.00
FIL3174	67982510	No	PTFE	1.0	1000	£966.00
FIL5024	67862502	No	PP	0.2	50	£130.00
FIL5026	67882502	No	PP	0.2	200	£427.00
FIL5028	67862504	No	Depth PP	0.45	50	£130.00
FIL5030	67882504	No	Depth PP	0.45	200	£424.00
FIL5032	67832510	No	GMF GD1*	1.0	100	£237.00
FIL5034	67832520	No	GMF GD2*	2.0	100	£235.00



\* Particle retention rating