

Sartorius Secura® Precision Balances

SARTORIUS

Reduce your risk with integrated and intelligent protection systems. Secura® gives you the security and peace of mind of knowing that you have done everything right. Besides providing highly accurate weighing results and operating convenience, Also features built-in protection systems for complete reliability and regulatory compliance – the safe and secure way. Real-time guidance prompts for levelling, automatic internal adjustment, and 100% traceable, clear documentation with sample and batch identifiers make your lab work more efficient.

- Integrated protection systems for minimising your risk of obtaining incorrect results
- Real-time alarm messaging and interactive user guidance for levelling
- Automatic detection and active monitoring to ensure you are within the acceptable USP operating range for reliable results
- Fully automatic internal adjustment function, isoCAL, for the highest possible accuracy at all times
- Cal Audit Trail for easy, traceable documentation of all calibration and adjustment data
- Safe and clear documentation by sample and batch identifiers for keeping perfect track of weighing data
- GLP/GMP compliant printouts
- Standard password security features for protection against tampering and unauthorised access



Code	Alt Ref	Capacity, g	Readability, g	Plate Ø, mm	Windshield	Dims, w x d x h, mm	Price
BAL7008	SECURA213-1S	210	0.001	120	Manual	216 x 360 x 320	£2540.00
BAL7012	SECURA313-1S	310	0.001	120	Manual	216 x 360 x 320	£2690.00
BAL7016	SECURA513-1S	510	0.001	120	Manual	216 x 360 x 320	£3270.00
BAL7020	SECURA1102-1S	1100	0.01	180	No	216 x 360 x 95	£2280.00
BAL7024	SECURA2102-1S	2100	0.01	180	No	216 x 360 x 95	£2480.00
BAL7028	SECURA3102-1S	3100	0.01	180	No	216 x 360 x 95	£2640.00
BAL7032	SECURA5102-1S	5100	0.01	180	No	216 x 360 x 95	£3410.00

Sartorius Quintix® Precision Balances

SARTORIUS

Quintix® makes laboratory work more enjoyable. The intuitive touch-screen in combination with the most precise weighing technology make Quintix® laboratory balances the only choice to achieve the best results. Benefit from improved workflows by easy-to-use application programs. If you need to generate GLP-compliant documents, you will save considerable time and effort with “plug and work” connectivity.

- Intuitive operation of the balance without having to read a manual thanks to the new user interface and colour touch screen
- Integrated application programs and intuitive user guidance for greater efficiency in your daily lab routine
- Highest accuracy and convenience ensured at all times by the isoCAL fully automatic internal adjustment feature



Code	Alt Ref	Capacity, g	Readability, g	Plate Ø, mm	Windshield	Dims, w x d x h, mm	Price
BAL7044	QUINTIX213-1S	210	0.001	120	Manual	216 x 360 x 320	£2040.00
BAL7048	QUINTIX313-1S	310	0.001	120	Manual	216 x 360 x 320	£2359.00
BAL7052	QUINTIX513-1S	510	0.001	120	Manual	216 x 360 x 320	£2822.00
BAL7056	QUINTIX612-1S	610	0.01	180	No	216 x 360 x 95	£1620.00
BAL7060	QUINTIX1102-1S	1100	0.01	180	No	216 x 360 x 95	£1830.00
BAL7064	QUINTIX2102-1S	2100	0.01	180	No	216 x 360 x 95	£2040.00
BAL7068	QUINTIX3102-1S	3100	0.01	180	No	216 x 360 x 95	£2305.00
BAL7072	QUINTIX5100-1S	5100	1.0	180	No	216 x 360 x 95	£1520.00
BAL7076	QUINTIX5101-1S	5100	0.1	180	No	216 x 360 x 95	£1870.00
BAL7080	QUINTIX5102-1S	5100	0.01	180	No	216 x 360 x 95	£2973.00