



Sartorius Practum® Precision Balances

SARTORIUS

Get the best value for your money, without compromises in precision and reliability. Rely on consistent readings and excellent repeatability. Practum® has everything you need for standard applications.

- Accurate, consistent readings with excellent repeatability based on Sartorius quality and technology designed in Germany
- High reliability and durability provided by mechanical overload protection
- Extra-large feet and up-front level indicator for easy levelling
- Supervisor lock for preventing unintentional changes to the setup settings
- Effortless weighing thanks to the ergonomic and sturdy draft shield (0.001g versions only)

Code	Alt Ref	Capacity, g	Readability, g	Plate Ø, mm	Windshield	Price
BAL7088	PRACTUM213-1S	210	0.001	120	Manual	£1020.00
BAL7090	PRACTUM313-1S	310	0.001	120	Manual	£1160.00
BAL7092	PRACTUM612-1S	610	0.01	180	No	£901.00
BAL7094	PRACTUM1102-1S	1100	0.01	180	No	£978.00
BAL7096	PRACTUM2102-1S	2100	0.01	180	No	£1160.00
BAL7100	PRACTUM5101-1S	5100	0.1	180	No	£978.00
BAL7098	PRACTUM5100-1S	5100	1.0	180	No	£738.00



Our service options:

- 🔧 Fault Find & Fix
- 📄 Service Contracts
- 🔧 On-Site Repairs
- 🔧 Pipette Servicing

Clarity in Servicing



Our Knowledge + Your Equipment = Service Tailored to You

Email: servicecentre@scientific-labs.com

Telephone: 0115 977 5075