

OHAUS Adventurer® Precision Balances



Striking the ideal balance between inventive features and functional, uncomplicated weighing capabilities, the OHAUS Adventurer incorporates all of the applications necessary for routine weighing and measurement activities. With a colour touchscreen, GLP/GMP compliance capabilities, and two USB ports.

- Wide viewing angle full-colour touch screen
- Dual USB ports with GLP/GMP tracking capabilities
- Available with or without Ohaus's AutoCal™ automatic calibration system
- 1mg models have a two piece, top-mounted draftshield with doors that open wide, providing ample room for sample placement
- 0.01g and 0.1g models feature large weighing pans



Code	Alt Ref	Model	Capacity, g	Readability, g	Calibration	Pan Dims, mm	Dims, w x d x h, mm	Price
BAL0015	30127960	AX223	220	0.001	Internal	∅ 130	354 x 230 x 340	£1455.00
BAL0016	30127961	AX223/E	220	0.001	External	∅ 130	354 x 230 x 340	£1200.00
BAL0017	30127962	AX423	420	0.001	Internal	∅ 130	354 x 230 x 340	£1605.00
BAL0018	30127963	AX423/E	420	0.001	External	∅ 130	354 x 230 x 340	£1350.00
BAL0019	30122614	AX523	520	0.001	Internal	∅ 130	354 x 230 x 340	£1915.00
BAL0020	30122615	AX523/E	520	0.001	External	∅ 130	354 x 230 x 340	£1660.00
BAL0054	30206079	AX422	420	0.01	Internal	171 x 195	354 x 230 x 100	£1170.00
BAL0056	30206080	AX422/E	420	0.01	External	172 x 195	354 x 230 x 100	£915.00
BAL0058	30206081	AX822	820	0.01	Internal	173 x 195	354 x 230 x 100	£1305.00
BAL0060	30206082	AX822/E	820	0.01	External	174 x 195	354 x 230 x 100	£1050.00
BAL0023	30127966	AX1502	1520	0.01	Internal	175 x 195	354 x 230 x 100	£1435.00
BAL0024	30127967	AX1502/E	1520	0.01	External	175 x 195	354 x 230 x 100	£1180.00
BAL0025	30127968	AX2202	2200	0.01	Internal	175 x 195	354 x 230 x 100	£1465.00
BAL0026	30127969	AX2202/E	2200	0.01	External	175 x 195	354 x 230 x 100	£1210.00
BAL0027	30122616	AX4202	4200	0.01	Internal	175 x 195	354 x 230 x 100	£1590.00
BAL0028	30122617	AX4202/E	4200	0.01	External	175 x 195	354 x 230 x 100	£1335.00
BAL0029	30122618	AX5202	5200	0.01	Internal	175 x 195	354 x 230 x 100	£1915.00
BAL0062	30206083	AX2201	2200	0.1	Internal	175 x 195	354 x 230 x 100	£1125.00
BAL0064	30206084	AX2201/E	2200	0.1	External	175 x 195	354 x 230 x 100	£870.00
BAL0066	30122619	AX8201	8200	0.1	Internal	175 x 195	354 x 230 x 100	£1415.00
BAL0068	30122620	AX8201/E	8200	0.1	External	175 x 195	354 x 230 x 100	£1160.00

Accessories

Code	Alt Ref	Description	Price
BAL5162	80251395	Auxiliary display for TA	£46.00
BAL6434	12106964	PC cable 9 pin for Ohaus balances	£71.00
BAL3954	80850000	Anti theft device	£48.00
BAL3959	83020221	In use display cover for Pioneer	£16.00

Can't find what you're looking for?

Let us know! We have over 300,000 products available, so if you can't find it our website will. Visit today.

www.scientificlabs.co.uk

