

## LTE IC Cooled Incubators Cont.

Code	Alt Ref	Description	Capacity, L
INC2174	IC/300/SS	IC cooled incubator stainless steel external finish with solid doors	176
INC2176	IC/400/SS	IC cooled incubator stainless steel external finish with solid doors	294
INC2178	IC/600/SS	IC cooled incubator stainless steel external finish with solid doors	516
INC2180	IC/200/SG	IC cooled incubator stainless steel external finish with double glazed doors	92
INC2182	IC/300/SG	IC cooled incubator stainless steel external finish with double glazed doors	176
INC2184	IC/400/SG	IC cooled incubator stainless steel external finish with double glazed doors	294
INC2186	IC/600/SG	IC cooled incubator stainless steel external finish with double glazed doors	516

## IC Cooled Incubator Accessories

Code	Alt Ref	Description
INC2188	IC/ACC/01	18mm access port (side panels only)
INC2190	IC/HLT/01	High/low temperature limiter with alarm
INC2192	IC/LIT/01	Interior illumination with 24hr or 7 day timer
INC2194	IC/SCK/01	Internal splash proof 13 amp socket
INC2196	IC/SHF/01	Extra shelf for 92L model
INC2198	IC/SHF/02	Extra shelf for 176L model
INC2216	IC/SHF/03	Extra shelf for 294L model
INC2218	IC/SHF/04	Extra shelf for 516L model
INC2220	IC/STK/01	Stacking kit
INC2222	IC/STK/02	Stacking kit with base unit castors

## LMS Cooled Incubators – SERIES 1A



The LMS Series 1A Cooled Incubators offer reliable and accurate temperature-controlled cabinets, operating between -10°C to +50°C. With 4 model sizes to choose from these versatile cabinets will suit most workloads ranging from 76 to 272L capacity.

- Digital set/display
- Fan-assisted circulation
- Hermetically-sealed refrigeration system
- Self-evaporating condensation drip tray
- Preset overheat protection
- For 220/240V, single phase supply
- White pre-coated aluminium interior

A full range of upgrades including fruit fly modification are also available to suit your requirements. Please enquire for details.



## Series 1A – Standard, Single Set point Control

Code	Alt Ref	Capacity, L	No. of Shelves	Ext. Dims, w x d x h, mm
INC5300	80	76	2	520 x 580 x 750
INC5302	120	120	3	520 x 580 x 1070
INC5304	201	201	4	520 x 580 x 1430
INC5306	280	272	5	520 x 580 x 1750