



Microscope Cavity Slides

76 x 26mm slides with a spherically shaped cavity polished into the upper surface.

Supplied in packs of 50

Code	Description	Cavity Dims, Ø x d, mm
MIC3450	Single cavity slide	13 x 0.35
MIC3452	Double cavity slide	16 x 0.5

Microscope Cover Slips

Thicknesses	
No.0	0.085 - 0.130mm
No.1	0.130 - 0.155mm
No.1½	0.155 - 0.185mm
No.2	0.190 - 0.250mm

SQUARE in Packs of 200

Code	Dims, w x l, mm	Thickness
MIC3100	18 x 18	No.0
MIC3102	20 x 20	No.0
MIC3104	22 x 22	No.0
MIC3110	18 x 18	No.1
MIC3112	20 x 20	No.1
MIC3114	22 x 22	No.1
MIC3116	24 x 24	No.1
MIC3120	18 x 18	No.1½
MIC3122	20 x 20	No.1½
MIC3124	22 x 22	No.1½
MIC3126	24 x 24	No.1½
MIC3130	18 x 18	No.2
MIC3132	20 x 20	No.2
MIC3134	22 x 22	No.2
MIC3136	24 x 24	No.2

RECTANGULAR in Packs of 100

Code	Dims, w x l, mm	Thickness
MIC3200	22 x 26	No.0
MIC3202	22 x 32	No.0
MIC3204	22 x 40	No.0
MIC3206	22 x 50	No.0
MIC3208	22 x 64	No.0
MIC3220	22 x 26	No.1
MIC3222	22 x 32	No.1
MIC3224	22 x 40	No.1
MIC3226	22 x 50	No.1
MIC3228	22 x 64	No.1
MIC3230	24 x 32	No.1
MIC3232	24 x 40	No.1
MIC3234	24 x 50	No.1
MIC3240	22 x 26	No.1½
MIC3242	22 x 32	No.1½
MIC3244	22 x 40	No.1½
MIC3246	22 x 50	No.1½
MIC3248	22 x 64	No.1½
MIC3250	24 x 32	No.1½
MIC3252	24 x 40	No.1½
MIC3254	24 x 50	No.1½
MIC3260	22 x 26	No.2
MIC3262	22 x 32	No.2
MIC3264	22 x 40	No.2
MIC3266	22 x 50	No.2
MIC3268	22 x 64	No.2

CIRCULAR in Packs of 100

Code	Ø, mm	Thickness
MIC3270	10	No.0
MIC3272	11	No.0
MIC3274	12	No.0
MIC3276	13	No.0
MIC3278	15	No.0
MIC3280	16	No.0
MIC3282	18	No.0
MIC3284	19	No.0
MIC3286	22	No.0
MIC3288	24	No.0
MIC3290	25	No.0
MIC3292	32	No.0
MIC3300	10	No.1
MIC3302	11	No.1
MIC3304	12	No.1
MIC3306	13	No.1
MIC3308	15	No.1
MIC3310	16	No.1
MIC3312	18	No.1
MIC3314	19	No.1
MIC3316	22	No.1
MIC3318	24	No.1
MIC3320	25	No.1
MIC3322	32	No.1
MIC3330	10	No.1½
MIC3332	11	No.1½
MIC3334	12	No.1½
MIC3336	13	No.1½
MIC3338	15	No.1½
MIC3340	16	No.1½
MIC3342	18	No.1½
MIC3344	19	No.1½
MIC3346	22	No.1½
MIC3348	24	No.1½
MIC3350	25	No.1½
MIC3352	32	No.1½
MIC3366	13	No.2
MIC3372	18	No.2
MIC3374	19	No.2
MIC3376	22	No.2
MIC3378	24	No.2
MIC3380	25	No.2
MIC3382	32	No.2

