

Brand Transferpette® S Variable Volume Pipettors



The piston-operated pipette Transferpette® S provides maximum versatility and optimum quality over the entire volume range.

- Large, central pipetting button and separate ejection function
- Autoclavable at 121 °C (2 bar), acc. DIN EN 285
- Volume-change protection – 4-position volume display
- Easy Calibration – readjustment without special tools
- Short stroke of only 12.5mm to reduce the risk of RSI (Repetitive Strain Injury)

Single Channel

Code	Alt Ref	Volume Range, µL	Compatible Pipette Tips, µL
PIP5100	704768	0.1 - 1	20
PIP5132	704769	0.1 - 2.5	20
PIP5102	704770	0.5 - 10	20
PIP5104	704772	2.0 - 20	200
PIP5134	704773	5 - 50	200
PIP5106	704774	10 - 100	200
PIP5108	704778	20 - 200	200
PIP8039	704780	100 - 1000	1000
PIP5112	704782	500 - 5000	5000
PIP5114	704784	1000 - 10,000	10,000



8 Channel

Code	Alt Ref	Volume Range, µL
PIP8042	703700	0.5 - 10
PIP8044	703706	5 - 50
PIP8046	703708	10 - 100
PIP8048	703710	20 - 200
PIP8050	703712	30 - 300



12 Channel

Code	Alt Ref	Volume Range, µL
PIP8052	703720	0.5 - 10
PIP8054	703726	5 - 50
PIP8056	703728	10 - 100
PIP8058	703730	20 - 200
PIP8060	703732	30 - 300



Starter Kits

Items Supplied: 3 x Transferpette® S Variable Volume Pipettors, 3 x TipBox (filled), 3 x shelf/rack mounts.

Code	Alt Ref	Description	Includes
PIP5140	704790	Micro	D-1, D-10, D-100
PIP5142	704791	Midi	D-20, D-200, D-1000
PIP5144	704792	Macro	D-1000, D-5000, D-10,000
PIP5145	704793	Standard	D-10, D-100, D-1000

