

Gilson PIPETMAN® G Manual Pipettes



PIPETMAN® G offers one of the lowest pipetting forces commercially available. With a redesigned piston assembly system and high quality piston grease, PIPETMAN® G guarantees hours of comfortable pipetting. It relies on the legendary accuracy, precision, and reliability of PIPETMAN®, with improved pipetting comfort. Made of PVDF and stainless steel, PIPETMAN® G offers premium durability. The new clip ejector allows for convenient access to serviceable parts.

- Premium pipetting comfort
- Reliable results
- Built to last and convenient
- New multichannel models available – 4 eight and twelve channel pipettes

PIPETMAN® G starter kit includes:

- 3 PIPETMAN® G pipettes: P20G, P200G, P1000G
- 2 racks DIAMOND tips: D200, D1000
- 3 SINGLE™ pipette holders
- 'Gilson Guide to Pipetting' handbook
- 3 plastic tip ejectors

PIPETMAN® G microvolume kit includes:


- 3 PIPETMAN® G pipettes: P2G, P10G, P100G (stainless steel ejectors)
- 2 racks DIAMOND tips: DL10, D200
- 3 SINGLE™ pipette holders
- 'Gilson Guide to Pipetting' handbook
- 3 plastic tip ejectors
- Adapter for P2G and P10G tip ejectors

PIPETMAN® G four pipette kit

- 4 PIPETMAN® G pipettes: P2G, P20G, P200G, P1000G
- 3 racks DIAMOND tips: D10L, D200, D1000
- 4 SINGLE™ pipette holders
- 'Gilson Guide to Pipetting' handbook

PIPETMAN® G large volume pipette kit

- 2 PIPETMAN® G pipettes: P5000G, P10mLG
- 2 bags of 10 filters



Code	Description	Channel	Volume Range, µL	Price
F144054M	P2G	Single	0.2 - 2.0	£209.00
F144055M	P10G	Single	0.5 - 10	£209.00
F144056M	P20G	Single	2 - 20	£209.00
F144057M	P100G	Single	10 - 100	£209.00
F144058M	P200G	Single	20 - 200	£209.00
F144059M	P1000G	Single	100 - 1000	£209.00
F144066	P5000G	Single	500 - 5000	£209.00
F144067	P10mLG	Single	1000 - 10000	£209.00
F144068	P8x10G	8	1 - 10	£552.00
F144069	P12x10G	12	1 - 10	£645.00
F144070	P8x20G	8	2 - 20	£552.00
F144071	P12x20G	12	2 - 20	£645.00
F144072	P8x200G	8	20 - 200	£552.00
F144073	P12x200G	12	20 - 200	£645.00
F144074	P8x300G	8	30 - 300	£408.00
F144075	P12x300G	12	30 - 300	£403.00
F167900	Starter kit	Single	—	£570.00
F167800	Microvolume kit	Single	—	£570.00
F167360	Four pipette kit	Single	—	£740.00
F167920	Large volume kit	Single	—	£381.00