

## B&S Polarimeter Tubes and Accessories

Manufactured to high quality standards conforming to ICUMSA recommendations and are compatible with most makes of polarimeter.

### Standard Glass – 8mm

| Code    | Alt Ref | Description  |
|---------|---------|--|
| POL1006 | 35-46   | Centre fill, glass, central arm: length 100mm            |
| POL1008 | 35-47   | Centre fill, glass, central arm: length 200mm            |
| POL1010 | 35-57   | Centre fill, glass, cup-shaped filler: length 100mm      |
| POL1012 | 35-58   | Centre fill, glass, cup-shaped filler: length 200mm      |
| POL1034 | 35-10   | Metal end, centre fill, glass, central arm: length 100mm |
| POL1036 | 35-11   | Metal end, centre fill, glass, central arm: length 200mm |

### Flow & Temperature Control – 8mm

ADP400 water-jacketed. Glass.

| Code    | Alt Ref | Description                                |
|---------|---------|--|
| POL1038 | 36-57   | Funnel flow-through tube: length 100mm     |
| POL1040 | 36-58   | Funnel flow-through tube: length 200mm     |
| POL1042 | 36-67   | Continuous flow-through tube: length 100mm |
| POL1044 | 36-68   | Continuous flow-through tube: length 200mm |
| POL1046 | 36-77   | Centre fill: length 100mm                  |
| POL1048 | 36-78   | Centre fill: length 200mm                  |

### Slotted Lids

Bellingham + Stanley polarimeter tubes are manufactured to high quality standards conforming to ICUMSA recommendations and are compatible with most makes of polarimeter. Tube ends are precision ground with windows made from specially selected low strain glass in order to achieve highest accuracy optical rotation measurement.

| Code    | Alt Ref | Description   |
|---------|---------|---|
| POL1050 | 37-009  | ADP400 slotted lid for 36-78 polarimeter tube - 200mm             |
| POL1052 | 37-010  | ADP400 slotted lid for 36-77 polarimeter tube - 100mm             |
| POL1054 | 37-011  | ADP400 slotted lid for 36-58/36-59/36-68 polarimeter tube - 200mm |
| POL1056 | 37-012  | ADP400 slotted lid for 36-57/36-67 polarimeter tube - 100mm       |

### Water Pumps

| Code    | Description             |
|---------|-------------------------|
| FIL1104 | Plastic, Edwards        |
| FIL1110 | Plastic, BRAND          |
| FIL1102 | Borosilicate glass, MBL |
| FIL1106 | Plastic, Kartell        |

