

## Sartorius SafetySpace® Filter Tips Continued

SARTORIUS

Code	Alt Ref	Volume Range, µL	Tip Length, mm	Pack	Price
PIP8904	790011F	0.1 - 10	32	10 x 96	£135.00
PIP8905	783201	0.1 - 10+	46	10 x 96	£135.00
PIP8907	790021F	0.5 - 20	51	10 x 96	£135.00
PIP8911	790101F	21 - 120	51	10 x 96	£135.00
PIP8913	790201F	5 - 200	52.5	10 x 96	£135.00
PIP7877	790301F	5 - 300	52.5	10 x 96	£135.00
PIP7883	791211F	50 - 1200+	90	10 x 96	£150.00

+ Extended tip

## Eppendorf epT.I.P.S. Pipette Tips

eppendorf

epT.I.P.S.® - the original Eppendorf "Totally Integrated Pipetting System" - are optimally coordinated to Eppendorf pipettes and, of course, meet EN ISO 8655 requirements. This has resulted in minimal attachment and ejection forces with the highest level of tightness. epT.I.P.S.® pipette tips can be used with pipettes from other manufacturers. Optimal wetting properties, high transparency and special, certified purity grades are the visible expression of our product and production philosophy.

## epT.I.P.S.® Racks

- Biopur® grade
- Guaranteed PCR clean, sterile, pyrogen free, ATP -free, bacterial/human DNA free, bacterial/human DNase/RNase free

Code	Volume Range, µL	Tip Length, mm	Pack	Price
E0030075226	0.1 - 20	40	5 x 96	£53.00
E0030075234	2 - 200	53	5 x 96	£48.70
E0030075242	20 - 300	55	5 x 96	£48.70
E0030075250	50 - 1000	71	5 x 96	£48.70
E0030075269	50 - 1250	76	5 x 96	£48.70
E0030075277	50 - 1250	103	5 x 96	£48.70
E0030075285	250 - 2500	115	5 x 48	£43.40
E0030075293	100 - 5000	120	5 x 24	£48.50
E0030075307	200 - 5000	175	5 x 24	£57.00
E0030075315	500 - 10000	165	5 x 24	£58.00



## epT.I.P.S.® Boxes 2.0

- Colour coded trays for simple identification of tips and matching Eppendorf pipette
- Autoclavable for future use

Code	Volume Range, µL	Tip Length, mm	Pack	Price
E0030076125	0.1 - 10	34	96	£33.50
E0030076133	0.1 - 20	40	96	£33.50
E0030076141	0.5 - 20	46	96	£29.70
E0030076150	2 - 200	53	96	£29.70
E0030076168	20 - 300	55	96	£29.70
E0030076176	50 - 1000	71	96	£29.70
E0030076184	50 - 1250	76	96	£33.50
E0030076192	50 - 1250	103	96	£34.00
E0030076206	500 - 2500	115	48	£39.40
E0030076214	100 - 5000	120	24	£35.20

