



Glisseal[®] Greases

Silicone-free grease for lubricating ground glass joints, stopcocks etc. Supplied in 60g tube.

- Excellent thermal stability
- Chemically stable against most acids, alkalis and gases
- Insoluble in water
- Can be removed with laboratory detergent
- Only a very thin layer is required for excellent lubrication

Code	Description	Price
VAC1020	Normal grease	£38.00
VAC1022	High vacuum grease	£37.70



Stopcock Grease

Silicone Grease M494 is water-repellent, non-melting, tasteless and odourless silicone grease used for electrical insulation purposes; it may also be used in contact with potable water supplies. Silicone Grease M494 is approved in the United Kingdom to criteria laid down in the Ministry of Defence Standard 68-69 Issue 1 (Grease XG 250; NATO Code S-736).

Code	For	Pack	Price
VAC1100	General use	50g	£9.33



Apiezon Greases

- Apiezon L is an ultra high vacuum grease, halogen free, working temperature 10 to 30°C
- Apiezon M is a high vacuum grease, halogen free, working temperature 10 to 30°C
- Apiezon T is a medium temperature vacuum grease, working temperature 10 to 120°C
- AP101 is an anti seize vacuum grease, working temperature -40 to 180°C

Code	Description	Pack	Price
VAC1000	Apiezon L	50g	£132.00
VAC1002	Apiezon M	100g	£114.00
VAC1006	Apiezon T	25g	£122.00
VAC1010	Apiezon AP101	50g	£64.00



EO Labs
Microbiological Culture Media

See pages 600–604 or visit www.scientificlabs.co.uk

