

MAUSER® Wiva Containers

- Fit for use in various medical environments, also for hygienic autoclaving and incineration processes
- Available in three sizes: 30, 50 and 60L
- Failure-proof locking system
- Supplied with solid lids
- UN approved for class 6.2 infectious waste, UN 3291

Code	Capacity, L	Dims, w x d x h, mm	Weight, g
SAF7360	30	335 x 400 x 318	1,360
SAF7362	50	335 x 404 x 539	1,650
SAF7364	60	335 x 400 x 640	1,900



Bio-Bin® Clinical Waste Bins

econix⁺

Bio-bin® is a new concept for non sharp infectious waste disposal. Designed to reduce your costs and CO₂ emissions when disposing of infectious and other waste materials, the paper base construction significantly reduces the amount of plastic and weight that is sent for final waste.

- Fold flat design saves storage space
- Reduced cost when compared to plastic bins
- Designs reflect DOH waste segregation guidelines
- 30L now available with a "twist and lock" temporary seal
- Loop version with longer base to accommodate longer instruments
- UN3291 approved container, no additional packaging required for transport

Code	Colour	Capacity, L	Pack
SAF7405	Yellow	1	Each
SAF7294	Yellow	2	Each
SAF7616	Yellow (Loop version)	2.5	100
SAF7434	Yellow	5	100
SAF7296	Yellow	6	Each
SAF7618	Yellow	30	10
SAF7403	Orange	1	Each
SAF7412	Orange	2	100
SAF7432	Orange	5	6
SAF7442	Orange	6	40
SAF7406	Purple/Yellow	1	100
SAF7416	Purple/ Yellow	2	100
SAF7436	Purple/Yellow	5	100
SAF7446	Purple/Yellow	6	40



Accessories

Accessories for Econix Bio-bin® clinical waste bins

Code	Description
SAF7450	Z Stand