

Bio-Bin® Clinical Waste Bins

econix

Bio-bin® is a new concept for non-sharp infectious waste disposal. Designed to reduce your costs and CO₂ emissions when disposing of infectious and other waste materials, the paper base construction significantly reduces the amount of plastic and weight that is sent for final waste.

- Fold flat design saves storage space
- Constructed of 96% paper - often cheaper to dispose of compared to plastic bins
- Eco-friendly
- Each empty bin weighs just grams, so doesn't add more cost to your disposal bill
- Designs reflect DOH waste segregation guidelines
- UN3291 approved container, no additional packaging required for transport
- Fully autoclavable



Colour	Waste Type	For
Yellow	Anatomical and infectious waste	Incineration
Orange	Anatomical and infectious waste	Alternative treatment
Purple	Cytotoxic and/or cytostatic waste	Incineration
Tiger stripe	Non infectious offensive hygiene waste	Deep landfill
Plain	General/non infectious	

Code	Alt Ref	Colour	Capacity, L	Pack	Price
SAF7405	EXBB1.YL	Yellow	1	Each	£1.30
SAF7294	EXBB2.YL	Yellow	2	Each	£1.50
SAF7616	EXBBL00P.YL	Yellow (loop version)	2.5	100	£113.00
SAF7434	EXBB5.YL	Yellow	5	100	£152.00
SAF7296	EXBB6.YL	Yellow	6	Each	£8.00
SAF7618	EXBB30.YL	Yellow	30	10	£54.00
SAF7424	EXBB30.YL	Yellow	30	10	£59.00
SAF7403	EXBB1.OR	Orange	1	Each	£1.20
SAF7412	EXBB2.OR	Orange	2	100	£106.00
SAF7456	EXBBL00P.OR	Orange (loop version)	2.5	100	£122.00
SAF7432	EXBB5.OR	Orange	5	6	£15.30
SAF7442	EXBB6.YL	Orange	6	40	£81.00
SAF7422	EXBB30.OR	Orange	30	100	£59.00
SAF7406	EXBB1.PL	Purple	1	100	£89.00
SAF7416	EXBB2.PL	Purple	2	100	£112.00
SAF7436	EXBB5.PL	Purple	5	100	£152.00
SAF7446	EXBB6.PL	Purple	6	40	£81.00
SAF7430	EXBB5.TS	Tiger stripe	5	100	£152.00
SAF7440	EXBB6.TS	Tiger stripe	6	40	£81.00
EXBBL00P0	EXBBL00P.0	Plain (loop version)	2.5	100	£119.00
EXBBPIP64	EXBB6.4	Plain	6	40	£81.00
EXBB304	EXBB30.0	Plain	30	10	£59.00

Accessories

Code	Alt Ref	Description	Price
SAF7450	BEXBBZS60	Z Stand	£17.90

