

Stuart Overhead Stirrers Continued




SS20 High Performance Overhead Stirrer

A powerful stirrer designed for demanding laboratory applications. Suitable for stirring liquids up to medium viscosity, including oils and microbiological media, up to 25L capacity. Features a hollow shaft arrangement to facilitate easy adjustment of paddle height. Advanced overload protection and sophisticated, responsive speed control.

- Used with volumes up to 25L
- Low to medium viscosity liquids
- Hollow shaft
- Keyless chuck


| Code | Alt Ref | Description | Price |
|---------|---------|------------------------------------|---------|
| STI2333 | SS20 | Overhead stirrer, high performance | £880.00 |



Stands

Extra strong retort stands with H-pattern base for stability and robust support rod. Both versions include bosshead. Heavy duty stand recommended for use with SS30 overhead stirrer.

| Code | Alt Ref | Description | Rod Dims, Ø x L, mm | Price |
|---------|---------|------------------------|---------------------|---------|
| STI2334 | SS10/1 | Stand, general purpose | 16 x 700 | £264.00 |
| STI2337 | SS10/2 | Stand, heavy duty | 25 x 850 | £688.00 |



Stirrer Paddles

Made from high grade stainless steel, there is a choice of five paddle heads and two lengths of 8mm diameter paddle rod. The heads screw on and off the rods, so they can be mixed and matched depending upon stirring requirements. This gives maximum versatility and value for money.

| Code | Alt Ref | Description | Price |
|---------|---------|--------------------------------|--------|
| STI2340 | SS10/5 | Paddle rod, 350mm | £23.10 |
| STI2342 | SS10/6 | Paddle rod, 550mm | £28.40 |
| STI2344 | SS10/10 | Small paddle head | £36.80 |
| STI2346 | SS10/11 | Large paddle head | £34.70 |
| STI2348 | SS10/12 | Large propeller head (2 blade) | £36.80 |
| STI2350 | SS10/13 | Small propeller head (4 blade) | £34.70 |
| STI2352 | SS10/14 | Anchor paddle head | £34.70 |