


## Hamilton Removable Needles

Removable needles give you the benefit of a replaceable needle for the occasional bend or clog, and the minimal dead volume of a cemented needle. The unique removable needle design installs the needle precisely at the zero line. Standard needle length is 51mm.



### Small Hub

Compatible with removable needle syringes that are 100µL and smaller.

Supplied in packs of 3

Code	Alt Ref	Gauge	Description	Price
<b>SYR7334</b>	7758-03	22S	Sharp	£52.00
<b>SYR7342</b>	7770-01	22S	Blunt	£66.00
<b>SYR7336</b>	7758-01	22	Sharp	£53.00
<b>SYR7330</b>	7758-02	26S	Sharp	£52.00
<b>SYR7338</b>	7768-01	26S	Blunt	£67.00



### Large Hub

Compatible with removable needle syringes from 250µL and larger.

Supplied in packs of 6

Code	Alt Ref	Gauge	Description	Price
<b>SYR7332</b>	7779-02	26S	Sharp	£68.00
<b>SYR7340</b>	7780-01	26S	Blunt	£68.00
<b>SYR7344</b>	7780-03	22S	Blunt	£53.00



## Hamilton Luer Lock Needles

A large selection of needles that are compatible with Luer Lock syringes. Luer Lock is an industry standard tapered termination utilized by most syringe manufacturers including medical Hypodermic syringes. Metal hub (nickel plated brass) needles can be used with TLL syringes and LT or TLL connectors. The needles are 304 stainless steel. Standard needle length is 51mm.

Supplied in packs of 6

Code	Alt Ref	Gauge	Description	Price
<b>SYR7346</b>	90039	26S	Sharp	£78.00
<b>SYR7350</b>	90038	22S	Sharp	£78.00
<b>SYR7354</b>	91039	26S	Blunt	£79.00
<b>SYR7358</b>	91038	22S	Blunt	£79.00