



Microscope Cavity Slides

C21-59

76 x 26mm slides with a spherically shaped cavity polished into the upper surface.

Supplied in packs of 50.

Code	Description	Cavity Dims, Ø x d, mm
MIC3450	Single cavity slide	16 x 0.2
MIC3452	Double cavity slide	16 x 0.5
MIC3454	Single cavity, economy	12 x 0.2

Microscope Cover Slips

C21-59

Thicknesses	
No.0	0.085 - 0.130mm
No.1	0.130 - 0.155mm
No.1½	0.155 - 0.185mm
No.2	0.190 - 0.250mm

SQUARE in Packs of 200

Code	Dims, w x l, mm	Thickness
MIC3100	18 x 18	No.0
MIC3102	20 x 20	No.0
MIC3104	22 x 22	No.0
MIC3110	18 x 18	No.1
MIC3112	20 x 20	No.1
MIC3114	22 x 22	No.1
MIC3116	24 x 24	No.1
MIC3120	18 x 18	No.1½
MIC3122	20 x 20	No.1½
MIC3124	22 x 22	No.1½
MIC3126	24 x 24	No.1½
MIC3130	18 x 18	No.2
MIC3132	20 x 20	No.2
MIC3134	22 x 22	No.2
MIC3136	24 x 24	No.2

RECTANGULAR in Packs of 100

Code	Dims, w x l, mm	Thickness
MIC3200	22 x 26	No.0
MIC3202	22 x 32	No.0
MIC3204	22 x 40	No.0
MIC3206	22 x 50	No.0
MIC3208	22 x 64	No.0
MIC3220	22 x 26	No.1
MIC3222	22 x 32	No.1
MIC3224	22 x 40	No.1
MIC3226	22 x 50	No.1
MIC3228	22 x 64	No.1
MIC3230	24 x 32	No.1
MIC3232	24 x 40	No.1
MIC3234	24 x 50	No.1
MIC3240	22 x 26	No.1½
MIC3242	22 x 32	No.1½
MIC3244	22 x 40	No.1½
MIC3246	22 x 50	No.1½
MIC3248	22 x 64	No.1½
MIC3250	24 x 32	No.1½
MIC3252	24 x 40	No.1½
MIC3254	24 x 50	No.1½
MIC3260	22 x 26	No.2
MIC3262	22 x 32	No.2
MIC3264	22 x 40	No.2
MIC3266	22 x 50	No.2
MIC3268	22 x 64	No.2

