

MILLEX® SYRINGE FILTERS YOU CAN TRUST.

**HUGE
SAVINGS!**



Millex® Syringe Filters

Millex® syringe filters set the standard for reliable small volume filtration for all your Chromatography or Life Science applications.

Contact us for great pricing if you favourite Millex® filter is not listed.

Code	Membrane	Ø, mm	Pore Size, µm	Sterility	Connections	Pack	List Price, £	Offer Price, £
SLCR033NK	Hydrophilic PTFE	33	0.45	N	Female Luer Lok/Male Luer Slip	1000	1800.00	1315.00
SLGL0250S	MCE	25	0.22	Y	Female Luer Lok/Male Luer Lok	50	196.00	173.00
SLHA033SS	MCE	25	0.45	Y	Female Luer Lok/Male Luer Slip	50	179.00	129.00
SLGSV255F	MCE	33	0.22	Y	Female Luer Lok/Male Luer Slip	50	241.00	183.00
SLHA02510	MCE	33	0.45	N	Female Luer Lok/Male Luer Slip	1000	1730.00	885.00
SLHN033NK	Nylon	33	0.45	N	Female Luer Lok/Male Luer Slip	1000	1490.00	1310.00
SLGP033RS	PES	33	0.22	Y	Female Luer Lok/Male Luer Slip	50	202.00	65.00
SLGP033RB	PES	33	0.22	Y	Female Luer Lok/Male Luer Slip	250	710.00	236.00
SLGPR33RB	PES	33	0.22	Y	Female Luer Lok/Male Luer Slip	250	520.00	143.00
SLFGN25VS	PTFE	25	0.2	Y	Female Luer Lok/Needle	25	146.00	79.00
SLFG025LS	PTFE	25	0.2	Y	Female Luer Lok/Male Luer Slip	50	237.00	166.00
SLFG05010	PTFE	50	0.2	N	Stepped Hose Barb/Stepped Hose Barb	10	117.00	96.00
SLGV013SL	PVDF	13	0.22	Y	Female Luer Lok/Male Luer Slip	100	348.00	295.00
SLHVX13NK	PVDF	13	0.45	N	Female Luer Lok/Male Luer Slip	1000	1590.00	942.00
SLGV033RS	PVDF	33	0.22	Y	Female Luer Lok/Male Luer Slip	50	203.00	139.00
SLGV033RB	PVDF	33	0.22	Y	Female Luer Lok/Male Luer Slip	250	760.00	504.00
SLHVM33RS	PVDF	33	0.45	Y	Female Luer Lok/Male Luer Slip	50	202.00	159.00
SLHVR33RB	PVDF	33	0.45	Y	Female Luer Lok/Male Luer Slip	250	520.00	142.00
SLHV033NK	PVDF	33	0.45	N	Female Luer Lok/Male Luer Slip	1000	1710.00	617.00



MILLEX® SYRINGE FILTERS YOU CAN TRUST.

Sterile Syringe Filters for Research Applications

Tissue culture media and additives / Buffers / DMSO / Biological solutions. Available in sizes:

Membrane	Pore size (µm)	Diameter (mm)
PVDF Durapore® Membrane Lowest binding membrane for protein rich solutions	0.1 Sterile filtration & Mycoplasma removal	33
	0.22 Sterile filtration	33
		25
		13
		4
	0.45 Clarification of sterile solutions	33
		13
		4
	5.0 Clarification of sterile solutions	25
PES Millipore Express® PLUS Fast flow and low binding for cell culture media preparation	0.22 Sterile filtration	33
		25
	0.45 Clarification of sterile solutions	33
MCE MF-Millipore™ Most referenced general purpose membrane	0.22 Sterile filtration	33
		25
	0.45 Clarification of sterile solutions	33
		33
	0.80 Clarification of sterile solutions	25

Sterile Syringe Filters for Venting and Gas Filtration

Sterile filtering gases / Venting sterile containers /
In-line vacuum pump protection / Transducer protection.
Available in sizes:

Membrane	Pore size (µm)	Diameter (mm)
PTFE Fluoropore™ Membrane Hydrophobic chemistry for gas filtration	0.20	25
	0.20	50
	0.45	50
	1.0	50

Non-Sterile Syringe Filters for Chromatography

HPLC, IC, GC / General particle removal / Industrial /
Environmental.
Available in sizes:

Membrane	Diameter (mm)	Pore size (µm)
PES Millipore Express® Membrane Fastest flow, high throughput	13	0.22
		0.45
	33	0.22
		0.45
PVDF Durapore® Membrane Low-protein binding	4	0.22
	13	0.45
		0.22
	33	0.45
		0.22
	33	0.45
Nylon Membrane Broad chemical compatibility	13	0.20
		0.45
	33	0.20
		0.45
PTFE Hydrophilic Millipore® LCR Membrane Lowest extractables and excellent solvent resistance	4	0.20
		0.45
	13	0.20
		0.45
	33	0.20
		0.45


**SCIENTIFIC
LABORATORY
SUPPLIES**

United Kingdom:
slsiinfo@scientific-labs.com

Northern Ireland:
NI.sales@scientific-labs.com

Republic of Ireland:
ROI.sales@scientific-labs.com