

*Offers are valid until 31st December 2019. Conditions apply.

SLS Select Reagent Reservoirs

- Sterile versions are certified RNase, DNase and pyrogen free
- Drip-free pour spouts on all four corners for spill-free pouring
- Centre 'V' slot for maximum sample recovery
- Easy read fill line graduations
- 25mL dual reservoirs are divided into 10 and 15mL divisions



Code	Description	Pack	List Price, £	Offer Price, £
B2B6132	25mL sterile - individually wrapped	100	70.00	51.10
B2B6134	25mL sterile - inner packs of 5	200	85.00	63.00
B2B6136	25mL non-sterile - bulk pack	100	30.00	22.40
B2B6138	50mL sterile - individually wrapped	100	72.00	52.50
B2B6140	50mL sterile - inner packs of 5	200	88.00	65.10
B2B6142	50mL non-sterile - bulk pack	100	32.00	23.10
B2B6144	100mL sterile - individually wrapped	100	85.00	62.30
B2B6146	100mL sterile - inner packs of 5	200	100.00	73.50
B2B6148	100mL non-sterile - bulk pack	100	33.00	24.50

SLS Lab Pro Low Retention Tips

Low retention polymer designed to increase pipetting accuracy by virtually eliminating sample retention. Deliver volumes within 0.1% of the indicated volume, versus 0.7% for standard polymer tips thereby providing more accurate results.


Racked, Sterile, Filtered



Code	Description	Pack	List Price, £	Offer Price, £
B2B5654	10µL	10 x 96	74.00	33.00
B2B5656	10XL	10 x 96	74.00	33.00
B2B5658	20µL	10 x 96	74.00	33.00
B2B5660	100µL	10 x 96	74.00	33.00
B2B5662	200µL	10 x 96	74.00	33.00
B2B5664	300µL	10 x 96	74.00	33.00
B2B5666	1000µL	8 x 96	74.00	33.00
B2B5668	1250µL	8 x 96	74.00	33.00


Liquid Handling Offers – Down to the last drop

*Offers are valid until 31st December 2019. Conditions apply.



Reload, Non-Sterile

- Reload up to 10 trays in 90 seconds
- Fully recyclable wafers
- RNase, DNase and pyrogen free
- Requires 57% less storage space than full racks
- Reuse existing racks
- Generates 90% less packaging waste



Code	Description	Pack	List	Offer Price, £
B2B5670	10µL	10 x 96	33.36	16.00
B2B5672	10XL	10 x 96	31.84	15.50
B2B5674	200µL	10 x 96	31.84	14.00
B2B5676	1000µL	10 x 96	34.56	17.50

SLS Lab Basics Multi Channel Pipettes

- Smooth and light tip loading and ejection
- Leak free sealing of tips
- Manifold can rotate 360° for easy right or left hand operation
- Easily removable & fully autoclavable manifold
- 3 year warranty



Code	Channel	Volume Range, µL	Compatible Lab Pro tips	List Price, £	Offer Price, £
B2B5696	8	0.5 - 10	10µL/10XL	505.00	158.00
B2B5698	8	5 - 50	20µL/100µL	505.00	158.00
B2B5700	8	10 - 100	100µL/200µL	505.00	158.00
B2B5702	8	20 - 200	200µL/300µL	505.00	158.00
B2B5704	8	30 - 300	200L/300µL	505.00	158.00

Liquid Handling Offers – Down to the last drop

*Offers are valid until 31st December 2019. Conditions apply.

SLS Lab Basics Single Channel Pipettes



- Ergonomic design ensuring lightweight & soft plunger movement
- Click-stop system for volume selection
- Colour coded for easy identification
- Calibration conforms to EN-ISO 8655 standards ensuring high accuracy and precision
- Calibration report enclosed with every pipette
- 3 year warranty



Code	Volume Range, μL	Increment, μL	Compatible Lab Pro tips	List Price, £	Offer Price, £
B2B5678	0.1 - 2.5	0.01	10 μL /10XL	115.00	46.00
B2B5680	0.5 - 10	0.1	10 μL /10XL	115.00	46.00
B2B5682	2 - 20	0.1	10XL	115.00	46.00
B2B5684	5 - 50	0.5	20 μL /100 μL	115.00	46.00
B2B5686	10 - 100	0.5	100 μL /200 μL	115.00	46.00
B2B5688	20 - 200	1	200 μL /300 μL	110.00	46.00
B2B5690	100 - 1000	5	1000 μL /1250 μL	115.00	46.00
B2B5692	500 - 5000	50	--	138.00	52.00
B2B5694	1000 - 10000	100	--	138.00	52.00

SLS Lab Basics Bottle Top Dispensers with Re-circulating Valve



Ergonomically designed for easy use and outstanding performance in addition to being extremely durable and reliable. Calibrated in an ISO17025 accredited facility.

- Fully autoclavable
- Easy re-calibration: a single tool to recalibrate
- Volume lock: smooth and accurate volume setting
- Comes with 4 additional adapters for regular bottle sizes



Code	Volume, mL	List Price, £	Offer Price, £
B2B5806	1 - 10	180.00	108.00
B2B5808	2.5 - 30	190.00	114.00
B2B5810	5 - 60	190.00	114.00



Liquid Handling Offers – Down to the last drop

*Offers are valid until 31st December 2019. Conditions apply.

SLS Lab Basics Pipette Filler



Lightweight, UV resistant material. It features adjustable valves and is microprocessor controlled for reliability and accuracy. The powerful pump offers rapid yet sensitive aspiration and dispensing. Filling and dispensing speed is easily set using the same hand by pressing the large mode button. The bright LED's will indicate which of the three modes has been selected, high, low or gravity mode.

- For 1-100mL glass or plastic pipettes
- 5 different coloured autoclavable nose pieces as standard
- Uses standard available 0.45µm hydrophobic filters
- Long life, easily changeable rechargeable battery
- Sensitive and adjustable aspiration and dispensing
- Complete with universal voltage and plug faced charger
- Comes with bench stand and magnetic/fixable wall mount



Code	Description	List Price, £	Offer Price, £
B2B5706	SLS Lab Basics pipette filler	200.00	100.00

SLS Lab Pro Safevac Vacuum Aspiration System



Safevac is designed for laboratory waste recovery and separation of liquid and solid. It can be widely used in cell culture, DNA extraction, microplate waste removal, and any other liquid separation or recovery.

- Sensitive level sensor to prevent liquid overflow during operation
- Waste bottle, tubing and handle are fully autoclavable
- Lid, tubing and handle are easy to disassemble, clean and autoclave
- Wide range of adapters to fit tube, microplate and disk plate for various tasks



Code	Vacuum Range, mbar	Aspirating Speed, mL/s	Dims, w x d x h, mm	Bottle Volume, L	List Price, £	Offer Price, £
B2B5770	0 - 600	1 - 15	180 x 240 x 450	4	750.00	488.00