



PIPETMAN M

The Electronic Pipette With Guaranteed Results

The electronic PIPETMAN M has been designed with new easy-to-use features and guaranteed accuracy and precision in both standard and repetitive pipetting modes. Lightweight and comfortable, PIPETMAN M is ideal for even the most intensive pipetting routines. Currently with 33% off all PIPETMAN M models.



**33%
OFF**

PIPETMAN M SINGLE CHANNELS			
MODEL / VOLUMES	PART NUMBERS	LIST PRICE EACH	33% OFF
P10M (0.5-10µL)	F81022PR17	£383.00	£256.61
P20M (2-20µL)	F81023PR17		
P100M (5-100µL)	F81013PR17		
P200M (20-200µL)	F81024PR17		
P300M (20-300µL)	F81014PR17		
P1200M (100-1200µL)	F81015PR17		
P5000M (500-5000µL)	F81016PR17		
P10mLM (1-10mL)	F81017PR17		
PIPETMAN M MULTICHANNELS			
P8 x 10M (0.5-10µL)	F81025PR17	£758.00	£507.86
P8 x 20M (1-20µL)	F81027PR17		
P8 x 100M (10-100µL)	F81018PR17		
P8 x 200M (20-200µL)	F81029PR17		
P8 x 300M (10-300µL)	F81031PR17		
P8 x 1200M (50-1200µL)	F81020PR17	£875.00	£586.25
P12 x 10M (0.5-10µL)	F81026PR17		
P12 x 20M (1-20µL)	F81028PR17		
P12 x 100M (10-100µL)	F81019PR17		
P12 x 200M (20-200µL)	F81030PR17		
P12 x 300M (10-300µL)	F81032PR17		
P12 x 1200M (50-1200µL)	F81021PR17		

- **CUSTOMISABLE & SECURE:**
Advanced drag and drop software for creating pipetting protocols easily and rapidly.
- **EASY, ERGONOMIC PIPETTING:**
Simple two button operation, ergonomic design with virtually zero pipetting forces. Scroll effortlessly between five pipetting modes: Pipet (standard mode), Reverse, Repetitive, Mix, Custom.

Not to be used in conjunction with any other offer or contract discount. Valid until 31/08/17 to UK mainland customers only.

How To Order:

Call 0115 982 1111 or 01236 431857

for Scotland or email gilson@scientificlabs.com