

Corning® Lambda™ Plus Pipettors

Trusted design for your lab

CORNING



Lambda Plus single-channel pipettor



Lambda Plus 8-channel and 12-channel pipettors

Corning Lambda Plus Pipettors

Corning Lambda Plus pipettors have been engineered to provide the highest levels of accuracy, reproducibility, and comfort. Available in single-, eight- and twelve-channel configurations, all Lambda Plus pipettors feature smooth plunger movement and very low pipetting forces. Volume can be set by turning either the thumb wheel or the plunger button, and a locking mechanism prevents accidental volume change.

The Lambda Plus pipettors have been designed to maintain their calibration, even after repeated use and autoclaving. A key is included with each pipettor for easy in-house calibration.

- ▶ Factory calibrated to strict ISO 8655 standards
- ▶ 3-year warranty
- ▶ Fully autoclavable and UV resistant
- ▶ Designed for universal fit with all common brands of pipet tips
- ▶ Adjustable tip ejector

Corning Lambda Plus Starter Kit

The Lambda Plus starter kit contains four Lambda Plus single-channel pipettors conveniently packaged with useful accessories (manual, calibration key, color-coded rings, tube of piston grease, and identification stickers). The kit includes:

- ▶ Lambda Plus single-channel pipettor, 0.5 to 10 μL (Cat. No. 4071)
- ▶ Lambda Plus single-channel pipettor, 2 to 20 μL (Cat. No. 4072)
- ▶ Lambda Plus single-channel pipettor, 20 to 200 μL (Cat. No. 4074)
- ▶ Lambda Plus single-channel pipettor, 100 to 1000 μL (Cat. No. 4075)
- ▶ Universal carousel rack for six single-channel pipettors (Cat. No. 4078)
- ▶ Corning DeckWorks™ pipet tips: Cat. No. 4135 (10 μL), Cat. No. 4138 (200 μL), and Cat. No. 4140 (1000 μL)



Lambda Plus starter kit

Ordering Information

Corning® Lambda™ Plus Single-channel, 8-channel, and 12-channel Pipettors

Cat. No.	Volume Range (µL)	Channels	Accuracy (%)	Precision (%)	Non-filtered Tips (µL)	Qty/Cs
4070	0.2 - 2.0	1	±12 to ±1.5	<6.0 to <0.7	10	1
4071	0.5 - 10	1	±4.0 to ±0.5	<4.0 to <0.4	10	1
4072	2 - 20	1	±4.0 to ±0.8	<3.0 to <0.4	200	1
4073	10 - 100	1	±1.6 to ±0.8	<0.8 to <0.3	200	1
4074	20 - 200	1	±1.2 to ±0.6	<0.6 to <0.3	200	1
4075	100 - 1000	1	±1.0 to ±0.6	<0.45 to <0.25	1,000	1
4080	1 - 10	8	±8.0 to ±2.0	<6.0 to <1.2	10	1
4081	5 - 50	8	±4.0 to ±1.6	<2.5 to <0.6	200	1
4082	20 - 200	8	±3.0 to ±1.0	<1.5 to <0.6	200/300	1
4083	50 - 300	8	±1.6 to ±1.0	<1.5 to <0.6	300	1
4084	1 - 10	12	±8.0 to ±2.0	<6.0 to <1.2	10	1
4085	5 - 50	12	±4.0 to ±1.6	<2.5 to <0.6	200	1
4086	20 - 200	12	±3.0 to ±1.0	<1.5 to <0.6	200/300	1
4087	50 - 300	12	±1.6 to ±1.0	<1.5 to <0.6	300	1

Corning Lambda Plus Starter Kit

Cat. No.	Description	Qty/Cs
4069	Lambda Plus starter kit Includes: 4 single-channel pipettors: 0.5 to 10 µL (Cat. No. 4071), 2 to 20 µL (Cat. No. 4072), 20 to 200 µL (Cat. No. 4074), 100 to 1000 µL (Cat. No. 4075); universal linear stand for four single-channel pipettors (Cat. No. 6066); Corning DeckWorks pipet tips: 10 µL (Cat. No. 4135), 200 µL (Cat. No. 4138), 1000 µL (Cat. No. 4140)	1 kit

Accessories

4078	Universal carousel rack for six single-channel pipettors	1
4079	Universal linear rack for six pipettors (six single-channel, maximum of four 8-channel or 12-channel pipettors)	1
4088	Pipettor stand for one single-channel, 8-channel, or 12-channel pipettor	1
6066	Universal linear stand for four pipettors, orange (four single-channel, maximum of two 8-channel or 12-channel pipettors)	1
6067	Universal linear stand for four pipettors, transparent (four single-channel, maximum of two 8-channel or 12-channel pipettors)	1
6068	Universal linear stand for six pipettors, transparent (six single-channel, maximum of three 8-channel or 12-channel pipettors)	1

Warranty/Disclaimer: Unless otherwise specified, all products are for research use only. Not intended for use in diagnostic or therapeutic procedures. Corning Life Sciences makes no claims regarding the performance of these products for clinical or diagnostic applications.

For additional product or technical information, visit www.corning.com/lifesciences or call 800.492.1110. Outside the United States, call +1.978.442.2200 or contact your local Corning sales office.

CORNING

Corning Incorporated
Life Sciences

836 North St.
Building 300, Suite 3401
Tewksbury, MA 01876
t 800.492.1110
t 978.442.2200
f 978.442.2476

www.corning.com/lifesciences

ASIA/PACIFIC
Australia/New Zealand
t 61 427286832

China
t 86 21 3338 4338
f 86 21 3338 4300

India
t 91 124 4604000
f 91 124 4604099

Japan
t 81 3-3586 1996
f 81 3-3586 1291

Korea
t 82 2-796-9500
f 82 2-796-9300

Singapore
t 65 6572-9740
f 65 6735-2913

Taiwan
t 886 2-2716-0338
f 886 2-2516-7500

EUROPE
CSEurope@corning.com

France
t 0800 916 882
f 0800 918 636

Germany
t 0800 101 1153
f 0800 101 2427

The Netherlands
t 020 655 79 28
f 020 659 76 73

United Kingdom
t 0800 376 8660
f 0800 279 1117

All Other European Countries
t +31 (0) 206 59 60 51
f +31 (0) 206 59 76 73

LATIN AMERICA
grupoLA@corning.com

Brasil
t 55 (11) 3089-7400
Mexico
t (52-81) 8158-8400