## SAFE ASPIRATION STATION & KIT



## **SAFE ASPIRATION OF LIQUIDS**

Motorized Vacuum Unit for Quick and Easy Aspiration of Liquid Waste

Safe and Contamination-Free

Adaptable

Convenient









## Safe and Contamination-Free

The Safe Aspiration Station is a stand-alone motorized vacuum unit, ideal for ELISA tests, cell culture, and liquid waste disposal. The closed unit is designed to prevent contamination hazards, and all parts in contact with liquids are fully autoclavable.

The station is also equipped with a level detection system to prevent liquid overflow, and an in-line hydrophobic filter that traps aerosols for safe disposal of contaminated or hazardous liquids.

## Adaptable

Different adapters provided with the Safe Aspiration Kit can be fitted to the hand operator. The station also ensures an adjustable aspiration level indicated by a board of LEDs.

The choice of the adapter and the easy setting of the vacuum level enable quick and accurate aspiration of liquids from a variety of vessels.



Associated with the Safe Aspiration Kit, the Safe Aspiration Station is easy to use on a bench or in laminar flow cabinets. The unit is always ready to use, simply press the button of the aspiration hand operator to aspirate.

The ergonomic hand operator can be stored on the bottle handle. The Safe Aspiration Kit also contains a stand with storing gaps for the adapters and the hand operator to help your workspace remain clean and organized.









SAFE ASPIRATION STATION	
Dimensions	53 x 18 x 32 cm
$(W \times D \times H)$	(21 x 7 x 13 in.)
Weight	3.4 kg
Aspiration speed	8 L/min (air)
Aspiration rate	17 mL/s (with a Pasteur aspiration pipette)
Vacuum range	-300 to -600 mbar (adjustable)
Noise	< 50 dB(A) at 1 meter
Voltage	100-240 VAC, 50/60 Hz
Warranty	12 months

ORDERING INFORMATION		
DESCRIPTION	PART NUMBER	
SAFE ASPIRATION STATION		
Safe Aspiration Station, EU Plug	F110744	
Safe Aspiration Station, US Plug	F110745	
Safe Aspiration Station, UK Plug	F110746	
Safe Aspiration Station, AU Plug	F110747	
Safe Aspiration Station, JP Plug	F110748	
Safe Aspiration Station, CN Plug	F110749	
SAFE ASPIRATION KIT		
Safe Aspiration Kit	F110750	

www.gilson.com/ContactUs

