

## #

Air-Displacement Pipette

PREMIUM PIPETTING COMFORT



## **Pipetting Comfort**

PIPETMAN® G offers one of the lowest pipetting forces commercially available. With a redesigned piston assembly system and high quality piston grease, PIPETMAN G guarantees hours of comfortable pipetting.

## Reliable Results

PIPETMAN G relies on the legendary accuracy, precision, and reliability of PIPETMAN with improved pipetting comfort.

## Built-to-Last and Convenient

Made of PVDF and stainless steel, PIPETMAN G offers premium durability. The PIPETMAN G clip ejector allows for convenient access to serviceable parts.





	PIPETMAN® G Single Channel				
	Model	PIPETMAN® DIAMOND Tips		Part Number	Volume
er-the n	P2G	D10 DL10	DF10ST DFL10ST	F144054M	0.2-2 μL
PIS	P10G	D10 DL10	DF10ST DFL10ST	F144055M	1-10 μL
1, <sup>20</sup> / <sub>1/2</sub>	P20G	D200	DF30ST	F144056M	2-20 μL
\$ 100 <sub>E</sub>	P100G	D200	DF100ST	F144057M	10-100 μL
\$ 200 <sub>4</sub>	P200G	D200 D300	DF200ST DF300ST	F144058M	20-200 μL
9 100p	P1000G	D1000 D1200	DF1000ST DF1200ST	F144059M	100-1000 μL
PSOO	P5000G*	D5000		F144066	500-5000 μL
€ <sup>10</sup> % Pt0mL	P10mLG*	D10mL		F144067	1000-10000 μL
	PIPETMAN® G Multichannel				
210 kg	P8x10 P12x10	D10 DL10	DF10ST DFL10ST	F144068 F144069	1-10 μL
V-204¢ P20	P8x20 P12x20	DL10 D200	DFL10ST DF30ST	F144070 F144071	2-20 μL
8 200 gr	P8x200 P12x200	D200 D300	DF200ST DF300ST	F144072 F144073	20-200 μL
P300	P8x300 P12x300	D200 D300	DF200ST DF300ST	F144074 F144075	30-300 μL

Models equipped with a tip ejector are only sold with a stainless steel tip ejector. Plastic tip ejector is available as an accessory. Refer to PIPETMAN G User's Guide LT801122.

Gilson maximum permissible errors are guaranteed only when PIPETMAN pipettes are used with the recommended PIPETMAN DIAMOND Tips.

PIPETMAN G specifications are downloadable from www.gilson.com





ISO 9001 LT800591/B | ©2019 Gilson S.A.S. All rights reserved. All trademarks are the property of Gilson, Inc. and its subsidiaries.





<sup>\*</sup>P5000G and P10mLG do not have tip ejectors.