

# AmpliPur™

EXPERT TIPS



EXCELLENT FIT FILTER TIPS FOR MOLECULAR BIOLOGY APPLICATIONS

| Excellent Fit

| Ideal for Molecular Biology Applications


| ISO 8655 Compliant





## Excellent Fit\*

AmpliPur Expert Tips are designed to fit various pipette brands for maximum convenience on the bench. Equipped with a filter, they also help prevent cross-contamination.



## Ideal for Molecular Biology Applications

These tips significantly contribute to molecular testing performance by ensuring sample integrity during the pipetting process.

## ISO 8655 Compliant

Manufactured according to the certified Gilson quality standard, AmpliPur Expert Tips comply with ISO 8655 requirements and are certified free of RNase and DNase, preserving your sample integrity.



\*Pipette compatibility needs to be verified by the user.

Technical Information	
Gilson Quality System	ISO 9001
Product Compliance	ISO 8655
Analysis Certificate	Free of RNase and DNase Hazardous Materials (2015/863/EU)
Materials	Tips: Polypropylene (PP)
	Filter: Polyethylene (PE)
	Rack: Polypropylene (PP)
	Heat shrink: Plastic Optical Fibre (POF)
Packaging	10 racks of 96 tips (960 tips)
Disposable	Tips: Yes
	Rack: Reusable
Warranty	3 years

AmpliPur™ Expert Tips		
Volume Range	Nominal Volume of Compatible Pipettes**	Part Number (Commercial Units)
0.5-10 µL	2 µL / 10 µL	F174101
1-20 µL	20 µL	F174201
10-200 µL	100 µL*** / 200 µL 300 µL	F174301
100-1000 µL	1000 µL	F174401
100-1200 µL	1000 µL / 1200 µL	F174501

\*\*Within ISO 8655 specifications.

\*\*\*Compatible with many of 100 µL pipettes, except Gilson P100 single channel.

AmpliPur Expert Tips are only available in select countries. Please contact your local Gilson representative for more information.



[gilson.com/contactus](http://gilson.com/contactus)

**ISO 9001** LT800607/A | ©2021 Gilson S.A.S.  
**CERTIFIED** All rights reserved. All trademarks are the property of Gilson, Inc. and its subsidiaries.

