

# pipetman®

M connected



CONNECTED VARIABLE VOLUME  
ELECTRONIC PIPETTES



## | Make Manual Processes Verifiable

With PIPETMAN® M Connected, open new possibilities in protocol execution assistance and data acquisition through the Bluetooth® connection.

## | On-Target Performance

Trust in your results with guaranteed performance in both standard and repetitive pipetting mode.

## | Ease-of-Use

Save and transfer up to ten custom protocols from PIPETMAN M Software into your pipette. Executing protocols has never been so easy.



Pipette Smarter.  
Report Faster.  
Confidently.

 **GILSON®**



### PIPETMAN M Connected Single Channel

| Model  | Part Number | Volume Range        |                       | PIPETMAN DIAMOND Tips |                       |
|--------|-------------|---------------------|-----------------------|-----------------------|-----------------------|
|        |             | STANDARD PIPET Mode | REPETITIVE PIPET Mode |                       |                       |
| P10M   | F81040      | 0.5 - 10 µL         | 0.5 - 10 µL           | D10<br>DL10           | DF10ST<br>DFL10ST     |
| P20M   | F81041      | 2 - 20 µL           | 2 - 20 µL             | D200                  | DF30ST                |
| P100M  | F81042      | 5 - 100 µL          | 5 - 100 µL            | D200                  | DF100ST               |
| P200M  | F81043      | 20 - 200 µL         | 5 - 200 µL            | D200<br>D300          | DF200ST<br>DF300ST    |
| P300M  | F81044      | 20 - 300 µL         | 10 - 300 µL           | D200*<br>D300         | DF200ST*<br>DF300ST   |
| P1200M | F81045      | 100 - 1200 µL       | 20 - 1200 µL          | D1000*<br>D1200       | DF1000ST*<br>DF1200ST |
| P5000M | F81046      | 500 - 5000 µL       | 100 - 5000 µL         | D5000                 |                       |
| P10mLM | F81047      | 1 - 10 mL           | 200 µL - 10 mL        | D10mL                 |                       |

### PIPETMAN M Connected Multichannel

|           |        |              |              |               |                                |
|-----------|--------|--------------|--------------|---------------|--------------------------------|
| P8x10M    | F81048 | 0.5 - 10 µL  | 0.5 - 10 µL  | D10<br>DL10   | DF10ST<br>DFL10ST              |
| P12x10M   | F81049 |              |              |               |                                |
| P8x20M    | F81050 | 1 - 20 µL    | 1 - 20 µL    | DL10<br>D200  | DFL10ST*<br>DF30ST             |
| P12x20M   | F81051 |              |              |               |                                |
| P8x100M   | F81052 | 10 - 100 µL  | 5 - 100 µL   | D200          | DF100ST                        |
| P12x100M  | F81053 | 20 - 200 µL  | 5 - 200 µL   | D200<br>D300  | DF100ST*<br>DF200ST<br>DF300ST |
| P8x200M   | F81054 |              |              |               |                                |
| P12x200M  | F81055 | 10 - 300 µL  | 10 - 300 µL  | D200*<br>D300 | DF200ST*<br>DF300ST            |
| P8x300M   | F81056 |              |              |               |                                |
| P12x300M  | F81057 | 50 - 1200 µL | 50 - 1200 µL | D1200         | DF1200ST                       |
| P8x1200M  | F81058 |              |              |               |                                |
| P12x1200M | F81059 |              |              |               |                                |

Wirelessly connect your electronic pipette to PIPETMAN M Connected Software and Gilson applications.

Keep Track of Pipettes and Protocols Used

| Create your custom protocol on a PC or an Android tablet.

| Export your protocol steps for recording, reusing and sharing.

| Never miss a service date. Your pipette stores this information in the software for you to recall.



Gilson maximum permissible errors are guaranteed only when PIPETMAN pipettes are used with the recommended PIPETMAN DIAMOND Tips. Visit [www.gilson.com](http://www.gilson.com).

Under these conditions, Gilson volumetric specifications in standard pipetting (PIPET mode) are guaranteed with a performance exceeding ISO 8655-2 recommendations for this mode.

\*These PIPETMAN DIAMOND Tips can be used with the indicated pipettes for the REPETITIVE mode but not until the maximal volume of the pipette in PIPET mode. Please refer to the volume range of your tips.

Recommended tips to get best results in REPETITIVE mode.

PIPETMAN M Connected specifications are downloadable from [www.gilson.com](http://www.gilson.com)

[www.gilson.com/contactus](http://www.gilson.com/contactus)

**ISO 9001** LT800598/A | ©2019 Gilson S.A.S.  
CERTIFIED All rights reserved. All trademarks are the property of Gilson, Inc. and its subsidiaries.

