

GelBond® PAG Film

Polyacrylamide Support Film

GelBond® PAG Film is a transparent, flexible polyester film designed to support polyacrylamide or MDE™ Gels. The acrylamide monomers covalently attach to the coating on the film during the polymerization reaction. Gels remain permanently attached to the film through electrophoresis and all subsequent fixing, staining, destaining, and drying procedures.


■ Benefits

- **Reliable** – Polyacrylamide gels retain their original dimensions during staining and after drying
 - **Durable** – Gels, particularly thin ones, are easier to handle during staining, destaining, and drying when supported
 - **Convenient** – Gel orientation can be recorded directly on the GelBond® PAG Film prior to casting
- NOTE: Polyester films will not transmit light of less than 310 nm, and will fluoresce at higher wavelengths.

■ Applications

- Drying and support of polyacrylamide gels

☀ 18°C to 26°C, protect from light



GelBond® PAG Support Film

Ordering Information – GelBond® Film Sheets and Rolls

Cat. No. NA	Cat. No. EU	Product Name	Product Description	Sheet Size (mm)	Chamber Compatibility
54711	54711	GelBond® PAG Support Film Sheets	For polyacrylamide gels, 0.2 mm thick	138 mm × 158 mm (50 sheets)	
54723	54723	GelBond® PAG Support Film Sheets	For polyacrylamide gels, 0.2 mm thick	160 mm × 180 mm (50 sheets)	Hofer® SE400, SE600, Bio-Rad® PROTEAN® II
54727	54727	GelBond® PAG Support Film Sheets	For polyacrylamide gels, 0.2 mm thick	124 mm × 258 mm (50 sheets)	GE Multiphor®
54729	54729	GelBond® PAG Support Film Sheets	For polyacrylamide gels, 0.2 mm thick	220 mm × 165 mm (50 sheets)	
54731	54731	GelBond® PAG Support Film Sheets	For polyacrylamide gels, 0.2 mm thick	199 mm × 264 mm (50 sheets)	
54733	54733	GelBond® PAG Support Film Sheets	For polyacrylamide gels, 0.2 mm thick	203 mm × 260 mm (50 sheets)	GE Multiphor® II
54746	54746	GelBond® PAG Support Film Sheets	For polyacrylamide gels, 0.2 mm thick	350 mm × 430 mm (10 sheets)	X-ray size

Custom-cut GelBond® PAG Support Film is available upon special request. Please inquire for pricing and availability.