

Distribution Partner of

Reagecon



SCIENTIFIC
LABORATORY
SUPPLIES

Turbidity Standards

Reagecon's Turbidity Standards for ratio and non-ratio instruments are composed of suspended polymer microspheres.

These Turbidity Standards remove the handling, stability and accuracy problems associated with traditional Formazin Turbidity Standards.

These products are extremely accurate ($\pm 1\%$), NIST traceable and US EPA approved and have guaranteed stability throughout their entire shelf life, even after opening the bottle.

Reagecon also offers a range of Stabilised Formazin Turbidity Standards which are marketed under the Turbiform™ brand. Details of these are available at www.reagecon.com



- Extensive range 0.0 – 4000 NTU
- Accuracy $\pm 1\%$
- NIST traceable
- Well defined spherical shape
- US EPA approved for drinking water turbidity
- Non-toxic and Non-carcinogenic
- Long stable shelf life at all concentrations
- No special handling or disposal requirements
- Ready to use, our range covers the full turbidity measurement range
- Certificates of Analysis & Safety Data Sheets available online.

For further information or enquiries, please contact: **Scientific Laboratory Supplies Ltd:**

United Kingdom: slsinfo@scientific-labs.com

Republic of Ireland: ROI.sales@scientific-labs.com

Northern Ireland: NI.sales@scientific-labs.com

International: slsexport@scientific-labs.com