

Colour (Color) Standards

Reagecon manufactures the full range of ASTM, Saybolt, Platinum-Cobalt (Hazen), Gardner, European Pharmacopeia and United States Pharmacopeia Colour Standards for use with ASTM, APHA, ACS, European and United States Pharmacopeia standard methods. The ASTM standard methods include D1500, D6045 and D1209. The products can be used to calibrate, control, qualify and validate colour measurement instruments.

- ASTM Colour Standard Sample A05- A7
- Saybolt Colour Standards S +30 to S -15
- Platinum-Cobalt Scale No. 0 - No. 1000
- Gardner Colour Standards GARD02 - GARD18
- European Pharmacopeia Standards (Opalescence, Primary and Standard Solutions)



- Produced in accordance with ASTM (D1500, D6045, D1209) APHA, ACS, EP & USP methods
- Consistency of product - Independent, Traceable, Certified
- For use as Calibration and/or Quality Control Standards
- Presented in high quality tamper evident bottles
- Various pack sizes available
- Ready to Use
- Customised Standards available
- Certificates of Analysis and Safety Data Sheets available online



For further information or enquiries, please contact: **Scientific Laboratory Supplies Ltd:**

United Kingdom: slsinfo@scientific-labs.com

Republic of Ireland: ROI.sales@scientific-labs.com

Northern Ireland: NI.sales@scientific-labs.com

International: slsexport@scientific-labs.com