

Superior Quality!

**volac** Glass

## Pasteur Pipettes



Soda Lime  
Glass

NEW:  
CERIUM SODA  
Glass

Borosilicate  
Glass

All **volac** Glass Pasteur Pipettes are precision made from carefully selected tubing on fully automatic machines ensuring total uniformity of production. Tips and tapers are centred and can be checked by rolling the pipettes on a flat surface. Length and taper of the tips are carefully controlled. Constrictions for plugging are neatly formed at a uniform distance from the fire polished end. Constrictions are heat formed which avoids weaknesses or contamination of the glass found when moulding or tooling processes are used (ISO7712).

**volac** Glass Pasteur Pipettes are made with a slightly heavier wall thickness to provide the best quality product allied to greater strength to give lower breakages during transit and reduced risk of breakage in use or when fitting PVC teats.

Glass Pasteur Pipettes are available either plain, pre-plugged with non-absorbent roving or plugged and sterilised in Gamma irradiation.

### Technical Details:

- available in standard sizes 150 mm and 230 mm, unplugged or plugged, sterilised or non-sterile
- made of soda lime glass, CERIUM SODA glass (NEW), borosilicate glass
- '50 dropper' pipettes also available
- not a medical device





# Glass Pasteur Pipettes, disposable

## Soda Lime Glass



Length	unplugged non-sterile soda lime glass	plugged non-sterile soda lime glass
	Art.-No.	Art.-No.
150 mm	D810	D810/PP
230 mm	D812	D812/PP

**Packing unit:** 1000 pcs (4 inner packs x 250 pcs)

## NEW: CERIUM SODA Glass

Length	plugged sterile CERIUM SODA glass
	Art.-No.
150 mm	D810/VA/CERIUM
230 mm	D812/VA/CERIUM

**Packing unit:** 250 pcs. (10 x 25)

## Borosilicate Glass

Length	unplugged non-sterile borosilicate glass	plugged non-sterile borosilicate glass
	Art.-No.	Art.-No.
150 mm	D810/BORO	D810/BORO/PP
230 mm	D812/BORO	D812/BORO/PP

**Packing unit:** 1000 pcs (4 inner packs x 250 pcs)

## Soda Lime Glass

Length	50-Dropper unplugged soda lime glass	50-Dropper plugged soda lime glass
	Art.-No.	Art.-No.
150 mm	D897	D897/PP

**Packing unit:** 1000 (4 inner packs x 250)

## Teat for Pasteur Pipettes

made of PVC, smooth for safety and durability

Art.-No.
D813

**Packing unit:** 100



## Poulten & Graf GmbH

Karl-Carstens-Str. 10,  
97877 Wertheim, Germany  
Tel.: +49 9342 9229-0  
www.poulten-graf.de  
sales@poulten-graf.com

All data are subject to change - Copyright:  
Poulten & Graf GmbH, GERMANY 2023/06

