



ART SoftFit~L pipette tips

Designed for optimal performance
with Rainin LTS style pipettes




ART SoftFit~L pipette tips

For Rainin LTS Pipettes


Thermo Scientific™ SoftFit~L™ pipette tips are designed to fit Rainin® LTS® style pipettes, providing the accuracy, precision and comfort you have come to expect when using your pipette.

- Achieve reliable results: Only SoftFit~L tips match the accuracy and precision of LTS-style tips
- Cylindrical collar and positive-stop design ensure a perfect fit, reducing attachment and ejection forces
- Certified to be free of RNase, DNase, DNA, ATP and endotoxin contamination
- 5 essential tip volumes: 20 µL, 200 µL, 300 µL, 1000 µL, and 1200 µL
- Color-coded inserts for easy matching of the tip and pipette volumes at a glance
- Tip inserts are labeled with lot and product part numbers for traceability
- Available with optional, low-retention technology

Cylindrical collar shape and positive-stop design for a perfect, ergonomic fit.







Color-coded inserts for easy tip volume identification



Minimize Waste and Increase Pipetting Accuracy

The hydrophobic surface of SoftFit~L low retention tips minimizes the amount of liquid remaining in the tip following dispensing, increasing accuracy and reproducibility. Low retention tips are ideal for use with viscous or detergent containing samples, nucleic acid and protein assays, or with expensive reagents.



Tip Type	Benefits	Recommended Use
SoftFit~L filtered tips 	<ul style="list-style-type: none"> • Blocks aerosols from contaminating the pipette and prevents cross-contamination of samples 	<ul style="list-style-type: none"> • Sensitive applications (cell culture, PCR, nucleic acid work) and samples • Volatile or corrosive samples
SoftFit~L non-filtered tips 	<ul style="list-style-type: none"> • Economic, high-quality tips 	<ul style="list-style-type: none"> • Non-critical applications and samples
SoftFit~L low retention tips 	<ul style="list-style-type: none"> • Hydrophobic surface maximizes sample delivery and conservation of expensive reagents • Improved pipetting accuracy • Reduced sample and reagent waste 	<ul style="list-style-type: none"> • Viscous samples • Samples containing detergents • Protein and DNA samples • Small volumes



Hinged Rack



Reload Insert



Reload Insert Tower

ART SoftFit-L Specifications and Ordering Information

Volume	Description	Packaging	Standard Tips	Low Retention Tips
Filtered Tips			Part Numbers	
SoftFit-L hinged racks: The versatile hinged rack provides one-handed access to tips and converts easily to a traditional lift-off lid. Color-coded inserts make it easy to identify tip volume at a glance.				
20µL	SoftFit-L 20µL sterile, filtered tips	96 tips/rack, 10 racks/pk, 5 pk/cs	2749-HR	2749-05-HR
200µL	SoftFit-L 200µL sterile, filtered tips	96 tips/rack, 10 racks/pk, 5 pk/cs	2769-HR	2769-05-HR
300µL	SoftFit-L 300µL sterile, filtered tips	96 tips/rack, 10 racks/pk, 5 pk/cs	2739-HR	2739-05-HR
1000µL	SoftFit-L 1000µL sterile, filtered tips	96 tips/rack, 8 racks/pk, 4 pk/cs	2779-HR	2779-05-HR
1200µL	SoftFit-L 1200µL sterile, filtered tips	96 tips/rack, 8 racks/pk, 4 pk/cs	2789-HR	2789-05-HR
SoftFit-L reload inserts: SoftFit-L reload inserts are an environmentally-friendly packaging option with a single tray of tips, for reuse of SoftFit-L hinged racks.***				
20µL	SoftFit-L 20µL sterile, filtered tip inserts	96 tips/insert, 10 inserts/pk, 5 pk/cs	2749-RI	2749-05-RI
200µL	SoftFit-L 200µL sterile, filtered tip inserts	96 tips/insert, 10 inserts/pk, 5 pk/cs	2769-RI	2769-05-RI
300µL	SoftFit-L 300µL sterile, filtered tip inserts	96 tips/insert, 10 inserts/pk, 5 pk/cs	2739-RI	2739-05-RI
1000µL	SoftFit-L 1000µL sterile, filtered tip inserts	96 tips/insert, 8 inserts/pk, 4 pk/cs	2779-RI	2779-05-RI
1200µL	SoftFit-L 1200µL sterile, filtered tip inserts	96 tips/insert, 8 inserts/pk, 4 pk/cs	2789-RI	2789-05-RI
***SoftFit-L reload inserts are for use with SoftFit-L empty racks				
Sterile, Non-Filtered Tips			Part Numbers	
SoftFit-L hinged racks: The versatile hinged rack provides one-handed access to tips and converts easily to a traditional lift-off lid. Color-coded inserts make it easy to identify tip volume at a glance.				
20µL	SoftFit-L 20µL sterile, non-filtered tips	96 tips/rack, 10 racks/pk, 5 pk/cs	3721-HR	3721-05-HR
200µL	SoftFit-L 200µL sterile, non-filtered tips	96 tips/rack, 10 racks/pk, 5 pk/cs	3751-HR	3751-05-HR
300µL	SoftFit-L 300µL sterile, non-filtered tips	96 tips/rack, 10 racks/pk, 5 pk/cs	3731-HR	3731-05-HR
1000µL	SoftFit-L 1000µL sterile, non-filtered tips	96 tips/rack, 8 racks/pk, 4 pk/cs	3781-HR	3781-05-HR
1200µL	SoftFit-L 1200µL sterile, non-filtered tips	96 tips/rack, 8 racks/pk, 4 pk/cs	3741-HR	—
SoftFit-L reload insert tower: an environmentally-friendly packaging option with reload trays of tips nested in a space saving tower to easily refill SoftFit-L racks*				
20µL	SoftFit-L 20µL sterile, non-filtered tip reload insert tower	96 tips/insert, 10 inserts/pk, 5 pk/cs	3721-RI	3721-05-RI
200µL	SoftFit-L 200µL sterile, non-filtered tip reload insert tower	96 tips/insert, 10 inserts/pk, 5 pk/cs	3751-RI	3751-05-RI
300µL	SoftFit-L 300µL sterile, non-filtered tip reload insert tower	96 tips/insert, 10 inserts/pk, 5 pk/cs	3731-RI	3731-05-RI
1000µL	SoftFit-L 1000µL sterile, non-filtered tip reload insert tower	96 tips/insert, 8 inserts/pk, 4 pk/cs	3781-RI	3781-05-RI
1200µL	SoftFit-L 1200µL sterile, non-filtered tip reload insert tower	96 tips/insert, 8 inserts/pk, 4 pk/cs	3741-RI	—

ART SoftFit-L Specifications and Ordering Information

Volume	Description	Packaging	Standard Tips	Low Retention Tips
Non-Sterile, Non-Filtered Tips			Part Numbers	
SoftFit-L hinged racks: The versatile hinged rack provides one-handed access to tips and converts easily to a traditional lift-off lid. Color-coded inserts make it easy to identify tip volume at a glance.				
20µL	SoftFit-L 20µL non-sterile, non-filtered tips	96 tips/rack, 10 racks/pk, 5 pk/cs	3722-HR	3722-05-HR
200µL	SoftFit-L 200µL non-sterile, non-filtered tips	96 tips/rack, 10 racks/pk, 5 pk/cs	3752-HR	3752-05-HR
300µL	SoftFit-L 300µL non-sterile, non-filtered tips	96 tips/rack, 10 racks/pk, 5 pk/cs	3732-HR	3732-05-HR
1000µL	SoftFit-L 1000µL non-sterile, non-filtered tips	96 tips/rack, 8 racks/pk, 4 pk/cs	3782-HR	3782-05-HR
1200µL	SoftFit-L 1200µL non-sterile, non-filtered tips	96 tips/rack, 8 racks/pk, 4 pk/cs	3742-HR	—
SoftFit-L in reload insert towers: SoftFit-L reload towers are an environmentally-friendly packaging option with reload trays of tips nested in a space saving tower to easily refill SoftFit-L racks*				
20µL	SoftFit-L 20µL non-sterile, non-filtered tip reload insert tower	96 tips/insert, 10 inserts/pk, 5 pk/cs	3722-RI	3722-05-RI
200µL	SoftFit-L 200µL non-sterile, non-filtered tip reload insert tower	96 tips/insert, 10 inserts/pk, 5 pk/cs	3752-RI	3752-05-RI
300µL	SoftFit-L 300µL non-sterile, non-filtered tip reload insert tower	96 tips/insert, 10 inserts/pk, 5 pk/cs	3732-RI	3732-05-RI
1000µL	SoftFit-L 1000µL non-sterile, non-filtered tip reload insert tower	96 tips/insert, 8 inserts/pk, 4 pk/cs	3782-RI	3782-05-RI
1200µL	SoftFit-L 1200µL non-sterile, non-filtered tip reload insert tower	96 tips/insert, 8 inserts/pk, 4 pk/cs	3742-RI	—

* 2 empty SoftFit-L racks per pack of towers

*** SoftFit-L reload inserts are for use with SoftFit-L empty racks

SoftFit-L Empty Racks

Description	Packaging	Part Number
SoftFit-L Empty Rack for 20µL, 200µL and 300µL tip inserts	10 racks, 5 pk/cs	3704-ER
SoftFit-L Empty Rack for 1000µL, 1200µL tip inserts	8 racks, 4 pk/cs	3706-ER

Find out more at thermofisher.com/pipettetips