

Stronghold+ (GM007*)

extended cuff · AQL 0.65 · enhanced chemical splash protection
silicone free · purple · medical grade · heavy weight · nitrile
unique Unigloves formulation · single use gloves

Physical Specifications:

Material: Nitrile (latex free)
Colour: Purple
Length: 300 mm
Weight: 7.3 g
Thickness: 0.12 mm
Finish: Powder free, chlorinated
Fit: Ambidextrous
Sterile: No

Order Codes and Size Information:

GM0071 | Extra small
GM0072 | Small
GM0073 | Medium
GM0074 | Large
GM0075 | Extra large

Packaging:

Each dispenser box contains 100 gloves
Each case contains 10 dispenser boxes

Classification:



AQL 0.65 (Pinholes) EN455 (Medical Grade)



Food Safe



EN ISO 374-1



EN ISO 374-5:
VIRUS

For more about the **Advanced** range, please **click on the image**



Features

- Free from silicone
- Tested to EN374
- Tested to ASTM D6978-05
- Extended cuff (300mm glove length)
- AQL 0.65
- Tested to EN1186
- Nitrile formulation conforms to the hand during use
- Beaded cuff
- Manufactured with textured palms and fingers
- Heavy-weight construction
- Touchscreen compatible

Benefits

- Silicone free formulation means no finger marks on metal and glass surfaces
- Resistant to a range of chemicals
- Resistant to a range of cytotoxic drugs when tested to ASTM D6978-05
- Extended cuff helps prevent liquids getting into the glove when working in wet conditions and provides additional protection from chemical splashes
- Ultra-low hole level; exceeds medical standards for pinholes
- Suitable for contact with all food types
- Nitrile formulation conforms to the hand during use for a more comfortable fit
- Beaded cuff provides extra strength during donning and helps prevent liquid roll back
- Textured palms and fingers provide improved grip, particularly when working in wet or oily conditions
- Heavy-weight construction is tough, durable and long lasting
- Suitable for use with mobile and tablet displays without removing gloves

Sectors

- Aerospace • Chemical manufacture and use • Dentistry
- Dental hygiene • Electronics • Engineering • Food processing
- Healthcare • Laboratory work • Manufacturing and assembly
- Pharmaceutical work • Printing



UNIGLOVES UK
Registered in England No: 4010200
VAT registration 791817787

3 Ambley Green, Gillingham Business Park,
Gillingham, ME8 0NJ, UK
+44 (0)800 049 6602
enquiries@unigloves.co.uk
unigloves.co.uk

UNIGLOVES GERMANY
Arzt- & Klinikbedarf
Handelsgesellschaft mbH

Camp-Spich-Straße 71
D-53842 Troisdorf-Spich,
+49 (0) 2241-9323-0
unigloves@unigloves.de
unigloves.de

