

<u>Spectrophotometer</u> Accessories



K

Jenway[®] Spectrophotometer Accessories

Jenway[®] spectrophotometers are very versatile and flexible instruments with an extensive range of accessories available which have been specifically designed to be easily interchangeable. The active accessories include the 8-cell automatic turret, 6 cell water heated automatic turret, sipper pump, Peltier and combined sipper/Peltier. Jenway also offer a range of passive accessories that can suit any application or budget.





Ordering information

| Model number | Description | For use with |
|--------------|--|--------------|
| 660201 | Sipper pump | 67 series |
| 660301 | Electronic Peltier for temperature controlled applications | 67 series |
| 660401 | Automatic 8 position cuvette holder | 67 series |
| 660501 | Water heated automatic six-cell changer | 67 series |
| 660701 | Combined Sipper/Peltier | 67 Series |
| 661201 | Single water heated 10x10mm cell holder | 67 Series |
| 685131 | Water heated 10x10mm single cell holder | 6850 |
| 685304 | Micro-cuvette holder | 6850 |
| 685204 | 10x10mm path length cuvette holder | 6850 |
| 685005 | 10 to 100mm cuvette holder | 6850 |
| 685401 | 8 position automatic cell changer | 6850 |



Ordering information

| Model number | Description | For use with |
|---------------|--|--------------------------|
| 630204 | 10x10mm cuvette holder | 63, 72, 73 and 67 series |
| 637071 | 16/24mm test tube holder | 63, 72, 73 and 67 series |
| 630005 | Variable path length cuvette holder | 63, 72, 73 and 67 series |
| 630304 | Micro-cuvette holder with reduced aperture | 63, 72, 73 and 67 series |
| 634001 | Manual 4 position cuvette holder | 63 Series |
| 648001 | Water heated single cuvette holder | 63 Series |
| 543001 | 40 column serial printer | 63 Series |
| 725201 | Heated cell 10x10 cuvette holder | 72 Series |
| 736201 | Water heated 10 x 10mm path length cuvette holder | 73 Series |
| 735401 | Automatic 8 position cuvette holder | 73 series |
| 735201 | Sipper pump | 73 series |
| 735301 | Electronic Peltier for temperature controlled applications | 73 Series |
| 735701 | Combined Sipper/Peltier | 73 series |
| 560101 | Internal Printer | 67 and 73 Series |
| SMP50/PRINTER | External Printer | 72 Series |

Cuvettes

Jenway offer a range of glass and quartz cuvettes as well as disposable plastic cuvettes. The glass and quartz cuvettes are available in a range of path lengths from 10 to 100mm. Jenway also offer a range of glass and quartz flow-through and micro-volume cuvettes. Glass cuvettes are suitable for measurements in the visible region whereas quartz cuvettes may be used in both the UV region and visible region.



| Part Code | Material | Description | Visible wavelength | UV/Visible wavelengths | Path length | Fill volume |
|-----------|----------|----------------------------------|-----------------------|---------------------------|-------------|---------------|
| 060 084 | Plastic | Pack of 100 plastic cuvettes | 1 | | 10mm | 2.4 to 4.5ml |
| 060 229 | Plastic | Pack of 500 plastic cuvettes | 1 | | 10mm | 2.4 to 4.5ml |
| 060 087 | Plastic | Pack of 100 plastic cuvettes | 1 | | 10mm | 1.5 to 3.0ml |
| 060 230 | Plastic | Pack of 100 plastic cuvettes | 1 | ✓ | 10mm | 2.5 to 4.5ml |
| 060 232 | Plastic | Pack of 500 plastic cuvettes | 1 | ✓ | 10mm | 2.5 to 4.5ml |
| 060 231 | Plastic | Pack of 100 UV plastic cuvettes | 1 | ✓ | 10mm | 1.5 to 3.0ml |
| 035 143 | Plastic | Pack of 100 plastic cuvettes, UV | 1 | ✓ | 10mm | 70µl to 1.5ml |
| 035 027 | Glass | 10mm path length cuvette | 1 | | 10mm | 2.0 to 3.5ml |
| 035 086 | Glass | 20mm path length cuvette | 1 | | 20mm | Up to 7.0ml |
| 035 029 | Glass | 40mm path length cuvette | 1 | | 40mm | Up to 14.0ml |
| 035 087 | Glass | 50mm path length cuvette | 1 | | 50mm | Up to 17.5ml |
| 035 079 | Glass | 100mm path length cuvette | 1 | | 100mm | Up to 35.0ml |
| 035 123 | Glass | Micro (500µl) cuvette | 1 | | 10mm | 500µl |
| 035 126 | Glass | Semi-micro (1ml) cuvette | 1 | | 10mm | 1ml |
| 035 045 | Glass | Flow-through cell | 1 | | 10mm | 80µl |
| 035 044 | Quartz | Flow-through cell | 1 | ✓ | 10mm | 1.8ml |
| 035 047 | Quartz | Flow-through cell | 1 | ✓ | 10mm | 80µl |
| 035 028 | Quartz | 10mm path length cuvette | 1 | ✓ | 10mm | 2.0 to 3.5ml |
| 035 056 | Quartz | 20mm path length cuvette | 1 | ✓ | 20mm | Up to 7.0ml |
| 035 030 | Quartz | 40mm path length cuvette | 1 | ✓ | 40mm | Up to 14.0ml |
| 035 281 | Quartz | 50mm path length cuvette | 1 | ✓ | 50mm | Up to 17.5ml |
| 035 282 | Quartz | 100mm path length cuvette | 1 | ✓ | 100mm | Up to 35.0ml |
| 035 124 | Quartz | UV ultra micro (50µl) cuvette | 1 | ✓ | 10mm | 50µl |
| 035 125 | Quartz | UV micro (500µl) cuvette | 1 | 1 | 10mm | 500µl |
| 035 127 | Quartz | UV semi-micro (1ml) cuvette | 1 | 1 | 10mm | 1ml |
| 035 138 | Quartz | UV ultra-micro (100µl) cuvette | 1 | 1 | 10mm | 100µl |
| | | | | | | |



A Bibby Scientific Brand

Bibby Scientific - UK (Group HQ) Beacon Road, Stone, Staffordshire, ST15 0SA, United Kingdom Tel: +44 (0)1785 812121 Fax: +44 (0)1785 813748 e-mail: sales@bibby-scientific.com

sales@bibby-scientific.com

www.jenway.com



Please scan the QR/Mobile Tag with your smartphone for more information

Find out more!!!