



# LABCOLD™

medical & scientific refrigeration



## ULTF301 -80°C Ultra Low Temperature Sparkfree Chest Freezer

Part Number	ULTF301
Capacity (Litres)	301
Style	Chest Freezer
Exterior Dimensions (HxWxD)	920 x 1300 x 650mm
Interior Dimensions (HxWxD)	650 x 1100 x 450mm
Weight	80kg
Operating Temperature	-65°C to -85°C
Interior Construction	PVC Coated Steel
Exterior Construction	Painted Steel
Defrost	Manual
Access Port	✓
Sparkfree Interior	✓
Open Door Alarm	✗
High Temperature Alarm	✓
Low Temperature Alarm	✓
Mains Failure Alarm	✓
Min/max Recording	Optional
Remote Alarm Contacts	✓
Digital Temperature Display	✓
Key Operated Lid Lock	✓
Internal Polystyrene Lids	✓
Castors	4
Vacuum Release Port	✓
Noise Level	≤51dba
Energy Consumption	7.9 kWh/24h
Eco Friendly Refrigerant	✓

Compact ultra low chest freezer featuring a spark free interior for laboratory safety.

Fitted with a lid lock as standard, this freezer also comes complete with a digital temperature display and castors for easy manoeuvrability.



Featuring an energy saving Danfoss compressor this freezer offers little temperature variation across the cabinet, making it ideal for the storage of samples which demand tight temperature tolerances.

Bespoke racking can be provided if required as can baskets at an extra cost. To increase the life of your freezer and to ensure you get optimum performance, preventative maintenance contracts are also available. Why not ask for details when you order?



01256 705 570

[www.labcold.com](http://www.labcold.com)

