

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: [www.sigmaldrich.com](http://www.sigmaldrich.com)Email USA: [techserv@sial.com](mailto:techserv@sial.com)Outside USA: [eurtechserv@sial.com](mailto:eurtechserv@sial.com)

## Product Specification

Product Name:

Zinc - granular, 20-30 mesh, ACS reagent, ≥99.8%

**Product Number:****243469**

CAS Number:

7440-66-6

MDL:

MFCD00011291

Formula:

Zn

Formula Weight:

65.39 g/mol

**Zn**

TEST	Specification
Appearance (Color)	Conforms to Requirements
Silver or Grey	
Appearance (Form)	Conforms to Requirements
Granules or Filings	
ICP Major Analysis	Confirmed
Confirms Zinc Component	
Complexometric EDTA	≥ 99.80 %
Iron (Fe)	≤ 0.01 %
Lead (Pb)	≤ 0.01 %
Suitability	Suitable
Suitable For Arsenic Determination	
Meets ACS Requirements	Current ACS Specification

Specification: PRD.2.ZQ5.10000041518

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at [Sigma-Aldrich.com](http://Sigma-Aldrich.com). For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.