

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
ProClin™ 300 - pkg of 2.0 L

Product Number: 48915-U

TEST	Specification
Appearance (Clarity) Free of Foreign Matter	Clear
Appearance (Color)	Colorless to Light Yellow
Appearance (Form)	Liquid
Color	≤ 3.0
VCS (Gardner) Scale	
Total Active Ingredient	3.0 - 3.6 %
RH-573 (MIT)	0.60 - 1.00 %
2-Methyl-4-Isothiazolin-3-One	
RH-651 (CMIT)	2.10 - 2.80 %
5-Chloro-2-Methyl-4-Isothiazolin-3-One	
Manufacture Date (MMYYYY)	
Registered Trademark	-----
ProClin is a trademark of the Dow Chemical Company (Dow) or an affiliated company of Dow	
Supplier Lot Number	Confirmed
Lot#	
Packaged From	Confirmed
Product# Lot#	

Specification: PRD.3.ZQ5.10000047157

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.