


## ANGLE ROTOR, 6-PLACE | 1556



### Rotor

max. RPM   max. RCF	9,000 min <sup>-1</sup>   10,595
max. capacity	6 x 94 ml
run-up   run-down, braked in sec	36   32
angle   max. noise level	35°   60 dB (A)
temperature in °C <sup>1)</sup>	0
<b>Cat. No.</b>	<b>1556</b>



Lid bioseal<sup>®</sup>, phenol-resistant

Cat. No.



**INCLUSIVE**



### Vessels

capacity in ml	1,5	2	15	50	50	85	7,5–8,5	9–10	10	8,5–10	15	50	10	30	50
Ø x L in mm	11 x 38	11 x 38	17 x 100	35 x 105	34 x 100	38 x 101	15 x 92	16 x 92	15 x 102	16 x 100	17 x 120	29 x 115	16 x 80	26 x 95	29 x 107
max. RCF <sup>2)</sup>	10,324	10,324	10,052	10,414	10,414	10,595	10,052	10,052	10,052	10,052	10,052	10,052	10,414	9,690	10,142
radius in mm	114	114	111	115	115	117	111	111	111	111	111	111	115	107	112
<b>Cat. No.</b>	<b>microliter tubes</b>			<b>tubes<sup>2)</sup></b>			<b>blood collection / urine vessels</b>				<b>tubes with screw cap</b>				



### Adapter

						-									
boring Ø x L in mm	11,4 x 39	11,4 x 39	17,5 x 78	35 x 89	35 x 89	-	17,5 x 78	17,5 x 78	17,5 x 78	17,5 x 78	17 x 106	29,8 x 97	16,5 x 72	26 x 85	29 x 92
vessels per rotor	24	24	6	6	6	6	6	6	6	6	6	6	12	6	6
<b>Cat. No.</b>	<b>1449</b>	<b>1449</b>	<b>1478</b>	<b>1463</b>	<b>1463</b>	<b>-</b>	<b>1478</b>	<b>1478</b>	<b>1478</b>	<b>1478</b>	<b>1466</b>	<b>1454</b>	<b>1477</b>	<b>1447</b>	<b>1446</b>

### Vessels

capacity in ml	85	94
Ø x L in mm	38 x 106	38 x 102
max. RCF <sup>2)</sup>	10,595	10,595
radius in mm	117	117
<b>Cat. No.</b>	<b>tubes with screw cap</b>	



### Adapter

boring Ø x L in mm	-	-
vessels per rotor	6	6
<b>Cat. No.</b>	<b>-</b>	<b>-</b>

All footnotes in the catalog can be found on the last page.