

# PATHOSORB ABSORBENT



COMPONENT  
PACKAGING



The Intelsius PathoSorb is a Super Absorbent Polymer (SAP) that ensures excellent absorption capacity to size ratio. Upon absorption during accidental spillage, PathoSorb acquires a 'gel like' consistency, locking in the spilled material for easy disposal. PathoSorb is lint free, reducing sterile contamination risk.

Category A and Category B transport regulations specify that absorbent material must be placed between primary and secondary containers in sufficient quantity to absorb the entire contents of the primary containers.

## FEATURES

- ▶ PathoSorb utilises SAP technology to provide effective absorbency whilst maximising available payload space.
- ▶ For multiple sample transportation, PathoSorb 4 & 6 Bay pouches provide a dual purpose; simultaneously protecting and separating samples.

For more information please contact Intelsius UK [cs@intelsius.com](mailto:cs@intelsius.com)



# PATHOSORB ABSORBENT



COMPONENT  
PACKAGING



## PATHOSORB ABSORBENT RANGE

Product Code	Product Name	Dimensions	Water Absorbent Capacity	Saline Absorbent Capacity
AB001	PathoSorb 50ml	150 x 75 mm 6 x 3 in	56 ml 1.9 fl oz	50 ml 1.7 fl oz
AB010	PathoSorb 100ml	160 x 140 mm 6.3 x 5.5 in	112 ml 3.9 fl oz	100 ml 3.5 fl oz
AB009	PathoSorb 200ml	225 x 175 mm 8.9 x 6.9 in	223 ml 7.8 fl oz	200 ml 7 fl oz
AB002	PathoSorb 300ml	305 x 220 mm 12 x 8.7 in	335 ml 11.7 fl oz	300 ml 10.5 fl oz
AB007	PathoSorb 600ml	470 x 295 mm 18.5 x 11.6 in	693 ml 24.3 fl oz	600 ml 21 fl oz
AB005	PathoSorb 4 Bay Pouch	225 x 125 mm 8.9 x 4.9 in	169 ml 5.9 fl oz	152 ml 5.3 fl oz
AB008	PathoSorb 6 Bay Pouch	325 x 125 mm 12.8 x 4.9 in	283 ml 8.3 fl oz	215 ml 7.5 fl oz

**UNITED KINGDOM:** DGP Intelsius Ltd  
1 Harrier Court, Airfield Business Park,  
Elvington, York, YO41 4EA,  
United Kingdom

t: +44 (0) 1904 607 390  
f: +44 (0) 1904 607 391  
cs@intelsius.com

**GERMANY:** DGP Intelsius GmbH  
Eulerweg 11  
64291 Darmstadt  
Germany

t: +49 6151 4936650  
f: +49 (0) 6151 4936651  
info@intelsius.de

**AMERICA:** DGP Intelsius LLC  
1399 Perry Rd, Suite 190,  
Plainfield, IN 46168,  
United States of America

t: +1 317 452 4006  
f: +1 317 290 4621  
americasorders@intelsius.com

**IRELAND:** DGP Intelsius Ireland Ltd.  
Unit 8. Coolcaslagh Industrial Estate,  
Coolcaslagh, Killarney  
Co. Kerry, V93 K006, Ireland

t: +353 (0)64 668 61 71  
salesireland@intelsius.com

**MALAYSIA:** DGP Intelsius Sdn Bhd  
Plot 11, Lorong Bemban 2, Estet  
Perindustrian Bemban, 31000  
Batu Gajah, Malaysia

t: +605 3666 399  
f: +605 3669 399  
compliance@intelsius.com

**INDIA:** DGP Intelsius Cold Chain Pvt. Ltd.  
Sun Srishti, C 10, Saki Vihar Road,  
Chandivali Powai. Mumbai - 400072  
Maharashtra - India

+91 22 49746536  
compliance@intelsius.com